3) Detection

★ In order to prevent contamination, the following procedure (1) or (2) is recommended for detection and data recording, which enables the detection in a closed tube.

(1) Real-time turbidity detection

Real-time turbidity detection can be carried out with the Loopamp Realtime Turbidimeter. For information about the equipment, visit the Eiken GENOME SITE (URL: http://loopamp.eiken.co.jp/e/). For the detailed operation of the equipment, refer to the instruction manual of the equipment.

(2) End-point turbidity detection

End-point turbidity detection can be conducted with the Loopamp End Point Turbidimeter. For more information, refer to the instruction manual. Notice; There is no relationship between the end-point turbidity values and the initial amount of the template.

(3) Visual fluorescence detection

Visual fluorescence detection can be achieved by using the Loopamp Fluorescent Detection Reagent (available for sale separately)⁶⁾. UV lamp (wavelength at 240~260nm or 350~370nm), protective goggle or glass board are required for fluorescence detection.

When UV lamp of wavelength around 230nm is used, negative sample may look like radiating fluorescence. Judgment should be done by comparing fluorescence of sample with that of positive and negative controls.

When output of UV lamp is too strong, negative control may look like radiating fluorescence. In such a case, take the UV lamp away from the reaction tube or change the angle of the reaction tube so that the difference between positive and negative controls becomes observable.

The incubation can be carried out in the Loopamp Turbidimeters (Realtime or End Point) or in the commercially available incubators (required temperature accuracy within $\pm 0.5^{\circ}\text{C}$ with hot bonnet.). Turbidity detection is also compatible with visual fluorescence detection using fluorescent reagents. For more information, refer to the package insert for Loopamp Fluorescent Detection Reagent.

(4) Electrophoresis

In order to avoid contamination, extra care should be taken when handling the amplification products during electrophoresis process.

- 0.5~2 μL of the reaction solution is applied for electrophoresis with the 2% agarose gel.
- The gel is stained with ethidium bromide (EtBr) or SYBR Green I.
- A ladder pattern can be observed in electrophoresis, since the amplified products consist of various size of inverted repeats of the target sequence on the same strands.

[Caution for use]

1. Handling the kit

1) This reagent kit should be stored at -20°C . To prevent the reagents from deterioration, only take out the necessary amount of reagents from the freezer before use (No decline was observed in the kit performance even after repeated freezing and thawing for 20 times in the quality control testing. But, in order to maintain the reagents performance, keep off unnecessary freezing and thawing).

2) Thaw the reagents at room temperature and keep them on ice for reagents preparation and later use. Before use, spin down the tubes to drop down the solution staying on the tube wall or on the cap, mix well the solution and spin down again. Notice that fierce mixing should be avoided as it can inactivate the *Bst* DNA Polymerase.

2. Caution for visual fluorescence detection

For the preparation of the sample solution, do not use the buffers containing chelating reagents such as TE buffer. If chelating reagent is added to the reaction solution, manganese ion binding with Calcein would be chelated so that the fluorescence light is released even no amplification take place. Besides, a sample containing a large amount of Ca, Zn or Fe ion might cause the false positive test result.

3. Handling reaction tubes

1) Only use the specified Loopamp Reaction Tube for turbidity or fluorescence detection. Other reaction tubes might have different optical transparency and can cause misjudgment.

2) Take full care when handling reaction tubes, as they are vulnerable to scratches or damages.

3) Check carefully to see if the reaction tubes have any crack or scratch before use. Crack or scratch on the tube might not only cause false judgment but also contaminate the equipment. If the tubes are broken inside the reaction block of the Loopamp Turbidimeter (Realtime or End Point), the reaction solution can spill inside the equipment and cause unrecoverable contamination and malfunctioning.

4) By comparing the solution volume in all tubes, check visually if proper amount of sample solution/master mix. has been dispensed into the reaction tube.

4. Caution for amplification reaction

Since bubbles in the solution will interfere the turbidity measurement and cause false judgment, try not to cause any bubble when mixing the master mix. and the sample solution. If bubbles are present, spin down to get rid of the bubbles.

5. Handling reaction tubes after use

- 1) The caps of the used reaction tubes should not be opened. Contamination of amplified products on other samples may not only cause false judgment of the test result but also pollute testing area. In this case, a correct test result may not be obtained unless pollution is completely removed.
- 2) Keep the cap of the used tube completely closed and dispose it, according to the relevant regulations and instructions, by incineration or after double bagging it with sealable vinyl bag. To prevent the amplified products from dispersing, do not conduct autoclave sterilization treatment for disposal.

[Caution for Handling]

- 1) LAMP reaction is very sensitive and even the slightest amount of amplified product tainted into the reaction might cause false result. Therefore, avoid this type of contamination and carry out the sample and reagent preparation in different clean benches. Conduct the detection of the amplification using turbidimeter (for end-point or real-time detection) with which the detection as well as the reaction can be accomplished in a tube with its cap closed. Also, to avoid the possible contamination from the amplified product when the product is taken out from the tube for the electrophoresis detection. Conduct electrophoresis in a room different from that for sample and reagent preparation.
- 2) When ultraviolet lamp is used for the fluorescence visual judgment, do not stare directly at the UV light. Since UV light is harmful to the eyes, even watching for a short period would irritate eyes and cause symptoms similar to conjunctivitis. Look at it through glass board or protective goggles.
- 3) When handling the sample, always abide by the biohazard counter measures⁸⁾.
- 4) Do not expose the Loopamp Reaction Tube, master mix. preparation tubes to UV light. A change in color or degeneration caused by ultraviolet lamp sometimes results in misjudgment.
- 5) This kit is designed for research use only.
- 6) If the operator does not have the experience or knowledge in the field of nucleic acid testing, there is a possibility of false judgment. Therefore, make sure that the kit is used under the supervision of the experienced and knowledgeable technicians.
- 7) Eiken Chemical Co., Ltd. does not bear any responsibility for false judgment or any consequential damage derived from the false judgment caused by non-capability problems such as operation error.
- 8) Use the kit before the expiration date, which is labeled on the outer box (Exp.Date).
- 9) The reagent tube is made of polypropylene and the main material for kit case is paper. The institution disposing the reagent tube and case should bear the responsibility and abide by the clinical waste disposal regulations, water pollution prevention law, and any other regulation related.

[Unit, Storage, Expiration, Code No.]

Product Name	Unit	Storage	Expiration	Code No.
Loopamp® DNA Amplification Kit	48tests	-20°C	LMP204	
	96tests		1 year	LMP205
	192tests			LMP206

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Loopamp®
LAMP (Loop-mediated Isothermal Amplification) 法

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Read this instruction carefully before use

For research use only

Rev. January 2008 (Ver.4)

Loopamp®
LAMP (Loop-mediated Isothermal Amplification) method

DNA Amplification Kit

【Characteristics】

LAMP (Loop-mediated Isothermal Amplification) method is a novel gene amplification method capturing the following characteristics: ① Only one enzyme is required and the amplification reaction proceeds under isothermal condition^{1), 2)}. ② It has extremely high specificity because of the use of 4 primers recognizing 6 distinct regions on the target., ③ It has high amplification efficiency and enables amplification within a shorter time., ④ It produces tremendous amount of amplified products which makes simple detection possible^{3), 4), 5), 6)}.

This kit is a simple reagent kit for the amplification of DNA by the LAMP method, and amplification of target DNA can be conducted by simply incubating specimen and all reagents provided at a constant temperature (60~65°C) for a fixed period of time (1 hr. for standard).

【Contents of the kit】

	48 tests	96 tests	192 tests
(1) 2 × Reaction Mix. ^{*2} (RM) ^{*1}	0.6 mL	1 tube	2 tubes
(2) <i>Bst</i> DNA Polymerase (<i>Bst</i> DNA Polymerase) ^{*1}	60 μL	1 tube
(3) Distilled Water (DW) ^{*1}	1.0 mL	1 tube	2 tubes
(4) Primer Mix.DNA ^{*4} (PM DNA) ^{*1}	25 μL	1 tube	—
(5) Positive Control DNA ^{*3 *4} (PC DNA) ^{*1}	0.1 mL	1 tube	—

*1: The notation on each reagent tube is shown in () .

	Tris-HCl (pH8.8)	KCl	MgSO ₄	(NH ₄) ₂ SO ₄	Tween20	Betaine	dNTPs
	40 mM	20 mM	16 mM	20 mM	0.2 %	1.6 M	2.8 mM each

*2: composition (2 ×)

Tris-HCl (pH8.8)	40 mM
KCl	20 mM
MgSO ₄	16 mM
(NH ₄) ₂ SO ₄	20 mM
Tween20	0.2 %
Betaine	1.6 M
dNTPs	2.8 mM each

*3: The control DNA is a plasmid DNA that is inserted with a Hind III fragment (6,557 bp) of lambda phage DNA.

*4: Primer Mix.DNA and Positive Control DNA are not included in the 96 tests and 192 tests package.

【Principle】

LAMP method is a novel isothermal nucleic acid amplification method using 4 kinds of primers and DNA polymerase with strand displacement activity. Among 4 kinds of primers, two of them are inner primers, whose 3' end and 5' end are respectively designed to be complementary to two different regions of the target sequence. When the primer-linked strand is replicated, resulting the 3' end structure to self-anneals to the complementary region in self-structure, it forms the loop structure at the end.

This stem-loop structure will allow the 3' end to initiate self-elongation and another inner primer to anneal to its loop region to synthesize new DNA strand with strand-displacement manner. Through the repetition of these processes, the amplification proceeds, and this enables amplification to continue under isothermal condition with only one kind of enzyme.

For further details of the LAMP method, refer to Eiken GENOME SITE (URL: <http://loopamp.eiken.co.jp/e/>).

【How to use】

1. Materials required but not provided

- Sterilized tubes for master mix. preparation (0.5 mL, 1.5 mL)
- Micropipettes (0.5~10 μL, 10~100 μL, 100~1,000 μL)
- Pipette tips with filter
- Loopamp Reaction Tube
- Aluminum rack for cooling tubes
- Ice (crushed ice) and ice box
- Centrifuge for micro-tubes
- Centrifuge for 8-connected tubes
- Vortex mixer

< For real-time turbidity detection >

- Loopamp Realtime Turbidimeter^{*5}

< For end-point turbidity detection >

- Loopamp End Point Turbidimeter^{*5}

< For visual fluorescence detection >