3) Detection

★ **In order to prevent contamination, the following procedure (1) or (2) is recommended for detection and data recording, which enables the detection in a closed tube.**

- (1) Real-time turbidity detection
Real-time turbidity detection can be carried out with the Loopamp Realtime Turbidimeter. For information about the equipment, visit the Eiken GENOME SITE（URL：http://loopamp.eiken.co.jp/e/）. For the detailed operation of the equipment, refer to the instruction manual of the equipment.
- (2) End-point turbidity detection
End-point turbidity detection can be conducted with the Loopamp EndPoint Turbidimeter. For more information, refer to the instruction manual. Notice；There is no relationship between the end-point turbidity values and the initial amount of the template.
- (3) Visual florescence detection
Visual florescence detection can be achieved by using the Loopamp Fluorescent Detection Reagent (available for sale separately)[®]. UV lamp (wavelength at 240～260nm or 350～370nm), protective goggle or glass board are required for florescence detection .
When UV lamp of wavelength around 230nm is used, negative sample may look like radiating fluorescence. Judgment should be done by comparing fluorescence of sample with that of positive and negative controls.
When output of UV lamp is too strong, negative control may look like radiating fluorescence. In such a case, take the UV lamp away from the reaction tube or change the angle of the reaction tube so that the difference between positive and negative controls becomes observable.
The incubation can be carried out in the Loopamp Turbidimeters (Realtime or End Point) or in the commercially available incubators (required temperature accuracy within±0.5℃：with hot bonnet.) .
Turbidity detection is also compatible with visual florescence detection using fluorescent reagents. For more information, refer to the package insert for Loopamp Fluorescent Detection Reagent.
- (4) Electrophoresis
In order to avoid contamination, extra care should be taken when handling the amplification products during electrophoresis process.
 - ・ 0.5～2μL of the reaction solution is applied for electrophoresis with the 2 % agarose gel.
 - ・ The gel is stained with ethidium bromide (EtBr) or SYBR Green I .
 - ・ A ladder pattern can be observed in electrophoresis, since the amplified products consist of various size of inverted repeats of the target sequence on the same strands.

【Caution for use】

1. Handling the kit

- 1) This reagent kit should be stored at —20℃. To prevent the reagents from deterioration, only take out the necessary amount of reagents from the freezer before use (No decline was observed in the kit performance even after repeated freezing and thawing for 20 times in the quality control testing. But, in order to maintain the reagents performance, keep off unnecessary freezing and thawing.) .
- 2) Thaw the reagents at room temperature and keep them on ice for reagents preparation and later use. Before use, spin down the tubes to drop down the solution staying on the tube wall or on the cap, mix well the solution and spin down again. Notice that fierce mixing should be avoided as it can inactivate the Enzyme Mix. (EM) .

2. Caution for visual florescence detection

For the preparation of the sample solution, do not use the buffers containing chelating reagents such as TE buffer. If chelating reagent is added to the reaction solution, manganese ion binding with Calcein would be chelated so that the fluorescence light is released even no amplification take place. Besides, a sample containing a large amount of Ca, Zn or Fe ion might cause the false positive test result.

3. Handling reaction tubes

- 1) Only use the specified Loopamp Reaction Tube for turbidity or florescence detection. Other reaction tubes might have different optical transparency and can cause misjudgment.
- 2) Take full care when handling reaction tubes, as they are vulnerable to scratches or damages.
- 3) Check carefully to see if the reaction tubes have any crack or scratch before use. Crack or scratch on the tube might not only cause false judgment but also contaminate the equipment. If the tubes are broken inside the reaction block of the Loopamp Turbidimeter (Realtime or End Point) , the reaction solution can spill inside the equipment and cause unrecoverable contamination and malfunctioning.
- 4) By comparing the solution volume in all tubes, check visually if proper amount of sample solution/master mix. has been dispensed into the reaction tube.

4. Caution for amplification reaction

Since bubbles in the solution will interfere the turbidity measurement and cause false judgment, try not to cause any bubble when mixing the master mix. and the sample solution. If bubbles are present, spin down to get rid of the bubbles.

5. Handling reaction tubes after use

- 1) The caps of the used reaction tubes should not be opened. Contamination of amplified products on other samples may not only cause false judgment of the test result but also pollute testing area. In this case, a correct test result may not be obtained unless pollution is completely removed.
- 2) Keep the cap of the used tube completely closed and dispose it, according to the relevant regulations and instructions, by incineration or after double bagging it with sealable vinyl bag. To prevent the amplified products from dispersing, do not conduct autoclave sterilization treatment for disposal.

【Caution for Handling】

1. LAMP reaction is very sensitive and even the slightest amount of amplified product tainted into the reaction might cause false result. Therefore, avoid this type of contamination and carry out the sample and reagent preparation in different clean benches. Conduct the detection of the amplification using turbidimeter (for end-point or real-time detection) with which the detection as well as the reaction can be accomplished in a tube with its cap closed. Also, to avoid the possible contamination from the amplified product when the product is taken out from the tube for the electrophoresis detection. Conduct electrophoresis in a room different from that for sample and reagent preparation.
2. When ultraviolet lamp is used for the fluorescence visual judgment, do not stare directly at the UV light. Since UV light is harmful to the eyes, even watching for a short period would irritate eyes and cause symptoms similar to conjunctivitis. Look at it through glass board or protective goggles.
3. When handling the sample, always abide by the biohazard counter measures[®].
4. Since RNA is vulnerable to RNase, handle with care. There is possible contamination of RNase from such as sample specimens (blood and urine, etc.), tissues, experimental tools, reagents, water, the operator’s saliva or perspiration. RNase is stable against heat and might maintain its activity even after autoclave sterilization. In order to avoid RNase contamination, the following points should be concerned.
 - ・ Designate exclusive operation area and tool for RNA treatment.
 - ・ Use sterilized disposable tubes and pipette tips.
 - ・ Use RNase free water (for example, DEPC treated water) .
 - ・ The operator must wear gloves and masks to prevent RNase contamination to the specimen from his or her saliva and perspiration.
5. Do not expose the Loopamp Reaction Tube, master mix, preparation tubes to UV light. A change in color or degeneration caused by ultraviolet lamp sometimes results in misjudgment.
6. This kit is designed for research use only.
7. If the operator does not have the experience or knowledge in the field of nucleic acid testing, there is a possibility of false judgment. Therefore, make sure that the kit is used under the supervision of the experienced and knowledgeable technicians.
8. Eiken Chemical Co., Ltd. does not bear any responsibility for false judgment or any consequential damage derived from the false judgment caused by non-capability problems such as operation error.
9. Use the kit before the expiration date, which is labeled on the outer box (Exp.Date) .
10. The reagent tube is made of polypropylene and the main material for kit case is paper. The institution disposing the reagent tube and case should bear the responsibility and abide by the clinical waste disposal regulations, water pollution prevention law, and any other regulation related.

【Unit, Storage, Expiration, Code No.】

Product Name	Unit	Storage	Expiration	Code No.
Loopamp [®] RNA Amplification Kit (RT-LAMP)	48tests	—20℃	1 year	LMP244
	96tests			LMP245
	192tests			LMP246

【References】

- 1) Notomi T. *et al.* : Nucleic Acids Research **28**, No.12, e63 (2000)
- 2) Nagamine K. *et al.* : Clin. Chem. **47**, No. **9**, 1742-1743 (2001)
- 3) Mori Y. *et al.* : Biochem. Biophys. Res. Commun. **289**, No.1,150-154 (2001)
- 4) Tomita N. *et al.* Abstract for The 73rd Annual Meeting of the Japanese Biochemical Society (2000)
- 5) Mori Y. *et al.* Abstract for the 23rd Annual Meeting of the Molecular Biology Society of Japan (2000)
- 6) Tomita N. *et al.* Abstract for the 26th Annual Meeting of the Molecular Biology Society of Japan (2003)
- 7) Nagamine K. *et al.* : Molecular and Cellular Probes **16**,No.3,223-229 (2002)
- 8) The guideline for the bio-safety and bio-hazard (by the Japanese Society for Bacteriology): Japanese Journal of Bacteriology **54**, No.3, 667-715 (1999)

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Loopamp[®] LAMP (Loop-mediated Isothermal Amplification) Ⅱ RNA 増幅試薬キット (RT-LAMP)

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RNA Amplification Kit (RT-LAMP)

【Characteristics】

LAMP (Loop-mediated Isothermal Amplification) method is a novel gene amplification method capturing the following characteristics: ① Only one enzyme is required and the amplification reaction proceeds under isothermal condition^{1), 2),}, ② It has extremely high specificity because of the use of 4 primers recognizing 6 distinct regions on the target., ③ It has high amplification efficiency and enables amplification within a shorter time., ④ It produces tremendous amount of amplified products which makes simple detection possible^{3), 4), 5), 6)}.

This simple-to-use amplification kit is able to amplify target RNA by the LAMP method (RT-LAMP). By employing the Enzyme Mix (containing reverse transcriptase and DNA polymerase), reverse transcription and DNA amplification can be accomplished at a constant temperature (60~65℃, usually 63℃) in a short time (1 hr. for standard) in one-step.

【Contents of the kit】

	48 tests	96 tests	192 tests
(1) 2 × Reaction Mix. * ² (RM) * ¹	0.6 mL	1 tube	2 tubes
(2) Enzyme Mix. * ³ (EM) * ¹	50 μL	1 tube	2 tubes
(3) Distilled Water (DW) * ¹	1.0 mL	1 tube	2 tubes
(4) Primer Mix.RNA * ⁴ (PM RNA) * ¹	30 μL	1 tube	—
(5) Positive Control RNA * ⁴ (PC RNA) * ¹	60 μL	1 tube	—

*¹ : The notation on each reagent tube is shown in ().

*² : composition (2 ×)

Tris-HCl (pH8.8)	40 mM
KCl	20 mM
MgSO ₄	16 mM
(NH ₄) ₂ SO ₄	20 mM
Tween20	0.2 %
Betaine	1.6 M
dNTPs	2.8 mM each

*³ : It is a mixture of *Bsf* DNA Polymerase and AMV reverse transcriptase.

*⁴ : Primer Mix.RNA and Positive Control RNA are not included in the 96 tests and 192 tests package.

【Principle】

LAMP method is a novel isothermal nucleic acid amplification method using 4 kinds of primers and DNA polymerase with strand displacement activity. Among 4 kinds of primers, two of them are inner primers, whose 3' end and 5' end are respectively designed to be complementary to two different regions of the target sequence. When the primer-linked strand is replicated, resulting the 3' end structure to self-anneals to the complementary region in self-structure, it forms the loop structure at the end.

This stem-loop structure will allow the 3' end to initiate self-elongation and another inner primer to anneal to its loop region to synthesize new DNA strand with strand-displacement manner. Through the repetition of these processes, the amplification proceeds, and this enables amplification to continue under isothermal condition with only one kind of enzyme. When RNA is used as template, by initially adding the reverse transcriptase, one-step process from cDNA generation to DNA amplification can be achieved.

For further details of the LAMP method, refer to Eiken GENOME SITE (URL : <http://loopamp.eiken.co.jp/e/>) .

【How to use】

1. Materials required but not provided

- Sterilized tubes for master mix. preparation (0.5 mL, 1.5 mL)
- Micropipettes (0.5~10 μL, 10~100 μL, 100~1,000 μL)
- Pipette tips with filter
- Loopamp Reaction Tube
- Aluminum rack for cooling tubes
- Ice (crushed ice) and ice box
- Centrifuge for micro-tubes
- Centrifuge for 8-connected tubes
- Vortex mixer

< For real-time turbidity detection >
○ Loopamp Realtime Turbidimeter *⁵

< For end-point turbidity detection >
○ Loopamp End Point Turbidimeter *⁵

< For visual fluorescence detection >
○ Loopamp Fluorescent Detection Reagent
○ Loopamp Realtime Turbidimeter, Loopamp EndPoint Turbidimeter, or incubator with hot bonnet (temperature accuracy within±0.5℃)

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- Heat block (for termination of the reaction) *⁵
- UV lamp (wavelength at 240~260nm or 350~370nm) *⁵
- Protective goggle or glass board

*⁵ : For the information about applicable instrument, reaction termination function and condition for UV irradiation, refer to Eiken GENOME SITE (URL : <http://loopamp.eiken.co.jp/e/>) .

2. Primer design

Appropriate primer design is essential for amplification by LAMP. Use exclusive primer designing software for the primer design. Refer to LAMP primer designing software, PrimerExplorer at the website: <http://primerexplorer.jp/e/>.

Using the highly purified primers, a rapid gene amplification can be performed and stable reproducibility of the amplification can be obtained. The first screening for appropriate LAMP primers might not necessarily require highly purified primers. However, after the primers are determined, it is recommended to use purified FIP and BIP through HPLC or better purification.

3. Reagents preparation

- Take out the reagents stored at -20℃, and thaw them at room temperature. Once the reagents are thawed, keep them on ice.
- Preparation of master mix. (Operate on ice) .
 - The following amount of the component is required for one reaction.

○ Sample reaction	
<Reagents>	<Amount>
2 × Reaction Mix. (RM)	12.5 μL
Primer : FIP	40 pmol
BIP	40 pmol
Loop-F * ⁶	20 pmol
Loop-B * ⁶	20 pmol
F3	5 pmol
B3	5 pmol
Enzyme Mix. (EM)	1.0 μL
Distilled Water (DW)	X μL (ad q.s.)
Total	20.0 μL/test

*⁶ : Loop primers are not necessarily required. However, the use of Loop primers shortens the amplification time by about 1/3⁷⁾.

○ Control reaction	
<Reagents>	<Amount>
2 × Reaction Mix. (RM)	12.5 μL
PrimerMix. RNA (PM RNA)	2.5 μL
Enzyme Mix. (EM)	1.0 μL
Distilled Water (DW)	4.0 μL
Total	20.0 μL/test

★ For visual fluorescence detection, also add 1 μL of the Loopamp Florescent Detection Reagent(available for sale separately) and maintain the total mixture amount at 20 μL.

★ When used in combination with Loopamp Primer Sets, follow the preparation instructions of each primer set to prepare the master mix.

- After dispensing, gently tap the tubes for a few times (hereinafter referred to as tapping) , or mix the solution by repeatedly inverting the tube, or mix by vortex mixer at about 1 second×3 times. After mixing well, centrifuge the tubes for a few seconds (hereinafter referred to as spin down) . And the mixture can be used as the master mix. for the reaction. Notice that too much mixing by the vortex mixer might inactivate the enzyme, and assure that vortexing is conducted at 1 second×3 times. The prepared master mix. should be used as soon as possible.

4. Operation procedure

- Mixing of master mix. and sample solution (Operate on ice)
 - Dispense 20 μL of the master mix. into each Loopamp Reaction Tube (available for sale separately) .
 - Add 5 μL of sample RNA to the master mix., and the volume of the solution should be 25 μL in total. For control reactions, use 5 μL of Distilled Water (DW) for negative control, and 5 μL of Positive Control RNA (PC RNA) for positive control. Mix the solution well by pipetting or tapping the tube with the cap closed and then spin down. Be careful not to cause air-bubbles when mixing.
- Amplification reaction
 - Set the reaction tubes in Loopamp Turbidimeter (Realtime or End Point) or the incubator with hot bonnet (temperature accuracy within ±0.5℃) , and incubate them at 60~65℃ for 30~60 minutes (The reaction condition is dependent upon the characteristics of the primer used, therefore examine the optimum condition beforehand) .
 - Inactivate the enzyme and terminate the reaction by incubating the mixture for 5 minutes at 80℃ or 2 minutes at 95℃.