- (3) After the enzyme inactivation, irradiate the reaction tube with UV lamp (wavelength at 240~260nm or 350~370nm) from the bottom and observe the fluorescence light from the side of the tube wearing protective goggles that can efficiently block ultraviolet light. When green fluorescence similar to positive control is observed, the sample is judged as positive. When no fluorescence like negative control is observed, the sample is judged as negative. The copy number does not necessarily correlate with the fluorescence intensity.

★ **For examples of detection cases using this reagent, refer to the products introduction page in Eiken GENOME SITE (URL ; <http://loopamp.eiken.co.jp/e/>).**

#### 【Caution for use】

##### 1. Handling the reagent

- This reagent kit should be stored at −20℃. To prevent the reagent from deterioration, only take it out from the freezer before use.
- Keep the reagent on ice for reagents preparation and later use. Before use, spin down the tubes to recollect the solution on the tube wall or on the cap, well mix the solution and spin down again.
- Read cautions for handling each reagent carefully before use.

##### 2. Caution for the sample preparation

For the preparation of the sample solution, do not use buffers containing the chelating reagents such as TE buffer. If chelating reagent is added to the reaction solution, manganese ion binding with Calcein will be chelated and the fluorescence light is released even when no amplification takes place. Also, a sample containing a large amount of Ca, Zn or Fe ion might cause false positive test result.

##### 3. Handling reaction tubes

- Only use the specified Loopamp Reaction Tube for fluorescence detection. Other reaction tubes might have different optical transparency and cause misjudgment.
- Take full care when handling reaction tubes since they are vulnerable to scratches or damages.
- Check carefully to see if reaction tubes have any crack or scratch before use. Crack or scratch on the tube might not only cause false judgment but also contaminate the equipment. If the tubes are broken inside the reaction block of the Loopamp Turbidimeter (Realtime or End Point), the reaction solution can spill inside the equipment and cause unrecoverable contamination and malfunctioning.
- By comparing the solution volume in all tubes, check visually if proper amount of sample solution/master mix. has been dispensed into the reaction tube.

##### 4. Caution for the detection

- Carry out the fluorescence visual detection after the amplification reaction and successive inactivation of the enzyme. If the inactivation procedure failed, it might cause false judgment. If necessary, save the detection result as picture format using digital camera.
- Never irradiate the reaction tube with UV light for long period or it might causes false judgment due to the increasing of the fluorescence background or the quenching of the fluorescence light.

##### 5. Handling reaction tubes after use

- The caps of the used reaction tubes should not be opened. Contamination of amplified products on other samples may not only cause false judgment of the test result but also pollute testing area. In this case, a correct test result may not be obtained unless pollution is completely removed.
- Keep the cap of the used tube completely closed and dispose it, according the relevant regulations or instructions, by incineration or after double bagging it with sealable vinyl bag. To prevent the amplified products from dispersing, do not conduct autoclave sterilization treatment for disposal.

#### 【Caution for Handling】

- LAMP reaction is very sensitive and even the slightest amount of amplified product tainted into the reaction might cause false result. Therefore, avoid this type of contamination and carry out the sample and reagent preparation in different clean benches.
- When ultraviolet lamp is used for the fluorescence visual judgment, do not stare directly at the UV light. Since UV light is harmful to the eyes, even watching for a short period would irritate eyes and cause symptoms similar to conjunctivitis. Look at it through glass board or protective goggles.
- When handling the sample, always abide by the biohazard counter measures<sup>9)</sup>.

- Do not expose the Loopamp Reaction Tube, master mix. preparation tubes to UV light. Change in color or degeneration in these materials caused by ultraviolet lamp sometimes results in misjudgment.
- This kit is designed for research use only.
- If the operator does not have the experience or knowledge in the field of nucleic acid testing, there is a possibility of false judgment. Therefore, make sure that the kit is used under the supervision of the experienced and knowledgeable technicians.
- Eiken Chemical Co., Ltd. does not bear any responsibility for false judgment or any consequential damage derived from the false judgment caused by non-capability problems such as operation error.
- Use the kit before the expiration date, which is labeled on the outer box (Exp.Date) .
- The reagent tube is made of polypropylene and the main material for kit case is paper. The institution disposing the reagent tube and case should bear the responsibility and abide by the clinical waste disposal regulations, water pollution prevention law, and any other regulation related.

#### 【Unit, Storage, Expiration, Code No.】

| Product Name                                       | Unit    | Storage | Expiration | Code No. |
|--|---------|---------|------------|----------|
| Loopamp <sup>®</sup> Fluorescent Detection Reagent | 96tests | −20℃    | 1 year     | LMP221   |

#### 【References】

- Notomi T. *et al.* : Nucleic Acids Research **28**, No.12, e 63 (2000)
- Nagamine K. *et al.* : Clin. Chem. **47**, No. **9**, 1742-1743 (2001)
- Mori Y. *et al.* : Biochem. Biophys. Res. Commun. **289**, No.1,150-154 (2001)
- Tomita N. *et al.* Abstract for The 73rd Annual Meeting of the Japanese Biochemical Society (2000)
- Mori Y. *et al.* Abstract for the 23rd Annual Meeting of the Molecular Biology Society of Japan (2000)
- Tomita N. *et al.* Abstract for the 26th Annual Meeting of the Molecular Biology Society of Japan (2003)
- Nagamine K. *et al.* : Molecular and Cellular Probes **16**,No.3,223-229 (2002)
- The guideline for the bio-safety and bio-hazard (by the Japanese Society for Bacteriology):Japanese Journal of Bacteriology **54**,No.3,667-715 (1999)

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Read this instruction carefully before use

For research use only

Rev. January 2008 (Ver.4)

# Fluorescent Detection Reagent

## 【Characteristics】

LAMP (Loop-mediated Isothermal Amplification) method is a novel gene amplification method capturing the following characteristics: ① Only one enzyme is required and the amplification reaction proceeds under isothermal condition<sup>1),2)</sup>., ② It has extremely high specificity because of the use of 4 primers recognizing 6 distinct regions on the target., ③ It has high amplification efficiency and enables amplification within a shorter time., ④ It produces tremendous amount of amplified products which makes simple detection possible<sup>3),4),5),6)</sup>.

This reagent allows the amplification reaction to be visually detected, using the ultraviolet (UV) lamp, by adding this reagent to the amplification reagent of Loopamp DNA Amplification Kit or Loopamp RNA Amplification Kit (RT-LAMP) before use.

## 【Contents of the kit】

|                                      |               |
|--------------------------------------|---------------|
| Fluorescent Detection Reagent (FD)*1 | 96 tests      |
| Fluorescent Detection Reagent (FD)*1 | 0.1 mL×1 tube |

\*1：The notation on reagent tube is shown in ( ) .

## 【Principle】

LAMP method can produce large amount of DNA amplification product as well as its by-product, pyrophosphate. Calcein included in the kit initially binds with manganese ion and thus remain quenched. When the amplification reaction proceeds, the manganese ion is deprived of Calcein by the generated pyrophosphate, which results in the emission of fluorescence. The free Calcein is apt to bind to magnesium ion in the reaction mixture so that it strengthens the fluorescence emission<sup>6)</sup>. This principle is applied in fluorescence detection of the LAMP reaction.

## 【How to use】

### 1. Materials required but not provided

- Loopamp DNA Amplification Kit or Loopamp RNA Amplification Kit (RT-LAMP)
- Sterilized tubes for master mix. preparation (0.5 mL, 1.5 mL)
- Micropipettes (0.5~10μL, 10~100μL, 100~1,000μL)
- Pipette tips with filter
- Loopamp Reaction Tube
- Aluminum rack for cooling tubes
- Ice (crushed ice) and ice box
- Centrifuge for micro-tubes
- Centrifuge for 8-connected tubes
- Vortex mixer
- Loopamp Realtime Turbidimeter, Loopamp End Point Turbidimeter\*2,3, or incubator with hot bonnet (temperature accuracy within±0.5℃)
- Heat block (for termination of the reaction)\*2
- UV lamp (wavelength at 240~260nm or 350~370nm)\*2
- Protective goggle or glass board

\*2：For the information about applicable instrument, reaction termination function and condition for UV irradiation, refer to Eiken GENOME SITE (URL：http://loopamp.eiken.co.jp/e/).

\*3：For the fluorescence visual detection, Loopamp Turbidimeter (Realtime or End Point) is used as an incubator. As for the detection, a UV irradiator is required.

## 2. Reagents preparation

- Take out the Fluorescent Detection Reagent stored at −20℃ (not frozen) and place the reagent directly on ice. Loopamp DNA Amplification Kit and Loopamp RNA Amplification Kit (RT-LAMP) should be thawed at room temperature before use and also kept them on ice.
- Preparation of master mix. (Operate on ice) .
  - Dispense the necessary amount of Fluorescent Detection Reagent (FD) and Loopamp DNA amplification reagents or Loopamp RNA amplification reagents into the master mix. preparation tube according to the table indicated below. Mix Fluorescent Detection Reagent (FD) , 2×Reaction Mix. (RM) , Primers and Distilled Water (DW) using vortex mixer, then add *Bst* DNA Polymerase or Enzyme Mix. (EM) .

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\*4：Loop primers are not necessarily required. However, the use of Loop primers shortens the amplification time by about 1/3<sup>7)</sup>.

- After dispensing, gently tap the tubes for a few times (hereinafter referred to as tapping) , or mix the solution by repeatedly inverting the tube, or mix by vortex mixer at about 1 second×3 times. After mixing well, centrifuge the tubes for a few seconds (hereinafter referred to as spin down) . And the mixture can be used as the master mix. for the reaction. Notice that too much mixing with the vortex mixer might inactivate the enzyme, and assure that vortexing is conducted at 1 second×3 times. The prepared master mix. should be used as soon as possible.

★ For the preparation of the sample solution, do not use the buffer containing the chelating reagent such as TE buffer. For details, refer to the 【Caution for use】 .

★ When used in combination with Loopamp Primer Sets, follow the preparation instructions of each primer set to prepare the master mix.

## 3. Operation procedure (Follow the procedure of Loopamp DNA Amplification Kit or Loopamp RNA Amplification Kit (RT-LAMP) )

1) Mixing of master mix. and sample solution (Operate on ice)

- Dispense necessary amount of master mix. into each Loopamp Reaction Tube according to the reagent to be used as in the following table. Then add sample DNA or RNA to the master mix. and the volume of the solution should be 25 μL in total.

|                                 | Master mix. | Sample solution |
|---------------------------------|-------------|-----------------|
| DNA Amplification Kit           | 23 μL       | 2 μL            |
| RNA Amplification Kit (RT-LAMP) | 20 μL       | 5 μL            |

- Mix the solution well by pipetting or tapping the tube with the cap closed and then spin down. Be careful not to cause air-bubbles when mixing.

2) Amplification reaction and its detection.

- Set the reaction tubes in Loopamp Turbidimeter(Realtime or End Point) or the incubator with hot bonnet (temperature accuracy within±0.5℃), and incubate them at 60~65℃ for 30~60 minutes (The reaction condition is dependent upon the characteristics of the primer used, therefore examine the optimum condition beforehand) .
- After the amplification reaction, terminate the reaction by inactivating the enzyme using the heat-block (Never fail to inactivate the enzyme) . The following different conditions are employed for the inactivation of the enzyme for DNA Amplification Kit and RNA Amplification Kit (RT-LAMP) .

|  |
|--|
| Loopamp DNA Amplification Kit ; at 80℃ for 2min or 95℃, 2min<br>Loopamp RNA Amplification Kit (RT-LAMP) ; at 80℃ for 5min or 95℃, 2min |
|--|