# **Quick Protocol**



## ScreenFect®mRNA Transfection Reagent

### **Package Contents**

Cat. No.	ScreenFect®mRNA	Dilution Buffer
S-5001-2	0.2 ml	10 ml
S-5001	1.0 ml	50 ml
S-5001-3	5 x 1.0 ml	5 x 50 ml

### **Storage Conditions**


Store ScreenFect® Reagents at 4°C. Do not freeze. For optimal long term activity, do not allow ScreenFect® Reagents to warm to room temperature each time it is used.

After several months of storage without using the reagent a slight precipitation might occur. If vortexed thoroughly, this has no influence on the performance of ScreenFect® Reagents.

#### **General Considerations**

For optimal results, amounts of ScreenFect®mRNA and nucleic acid (NA) need to be optimized for each cell type and each NA used. An optimization protocol is provided in our product manual which can be downloaded from the Incella homepage. We strongly recommend the One-Step transfection method for all of our products. For transfection of adherent cells, remove the used medium and mix fresh medium with the transfection complexes. Then add the mix to the cells.

## ScreenFect® Protocol: One-Step Transfection



# ScreenFect® Products

ScreenFect®A	Multipurpose reagent (most suitable for pDNA transfection, suitable for RNA applications) with ver low cytotoxicity.	
ScreenFect®A-plus	Multipurpose reagent with optimized formulation requiring less reagent per transfection.	
ScreenFect®siRNA	Specialized reagent for best performance in siRNA delivery.	
ScreenFect®mRNA	Optimized reagent for the delivery of mRNA.	

For additional information regarding ScreenFect®mRNA and other ScreenFect® Products, visit our homepage (www.screenfect.com) and view our product pages and instruction manuals.

# **Quick Protocol**



## Protocol for mRNA Transfection

	Component	Procedure for one well (96-well-plate)	96-well	24-well	6-well		
	Reagent Dilution	Dilute 0.25 µl of ScreenFect®mRNA in Dilution Buffer to a final volume of 10 µl and mix thoroughly.	0.25 μl reagent 10 μl dilution	1.5 µl reagent 40 µl dilution	5 μl reagent 120 μl dilution		
		Important: Vortex the reagent once per day of use. Add ScreenFect®mRNA reagent directly into supplied buffer with rapid pipette mixing or vortexing.					
2	mRNA Dilution	Dilute a total of 75 ng mRNA in Dilution Buffer to a final volume of 10 $\mu l. $	75 ng 10 μl dilution	300 ng 40 µl dilution	1000 ng 120 µl dilution		
		Tip: The amount of mRNA depends on the particular protein it encodes or function it has in the cell.					
3	Complex formation	Combine the dilutions of mRNA and ScreenFect®mRNA and mix immediately using 10 rapid pipette strokes.  Leave for 20 min at room temperature for complex formation.	20 µl complexes	80 µl complexes	240 μl complexes		
		Important: Do not vortex!					
4	Cell preparation & transfection	Add 80 $\mu l$ freshly detached and resuspended cells to the complexes and mix with pipette.	Add 80 µl cell suspension	Add 420 µl cell suspension	Add 1250 µl cell suspension		
		Tip: The time-saving reverse cell transfection method may not be suited for all cell types. To transfect adherent cells, first remove and discard medium from cells, then add 80 µl fresh culture medium to transfection complexes, mix with pipette and immediately apply to cells.					
5	Cell plating	Transfer the cells and complexes to one well of a 96-well plate.	Transfer cells with complexes to plate	Transfer cells with complexes to plate	Transfer cells with complexes to plate		

**Note:** This protocol is a guideline. Values are suitable for easy for transfect cell lines. This protocol does not replace optimization experiments. View our manual for instructions. Serum does not affect the performance of ScreenFect®mRNA but we recommend avoiding antibiotics.

Cells must be mycoplasma free, in exponential growth phase and have even plating density across the entire surface area.