



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

AMRESKO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Agarose, For DNA Separations

Code: 0710

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	9012-36-6		
Appearance:		Storage:	Room Temperature
Fine white powder		Grade:	BIOTECHNOLOGY GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
DNase	NONE	PASS
EEO(-Mr)	<= <0.13	PASS
Endonuclease/ligase inhibitory factors	NONE	PASS
Gel Strength (1.5%)	>= 1200 g/cm2	PASS
Gelling Range (1.5%)	36 - 39 Degrees C	PASS
Melting Range (1.5%)	87 - 89 Degrees C	PASS
Moisture	<= <10 %	PASS
Protease	NONE	PASS
RNase	NONE	PASS
Sulfate	<= 0.15 %	PASS

Spec Set: 0710BIO

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed:

05/18/2016

Page 1 of 1