

Corporate Headquarters AMRESCO, LLC 6681 Cochran Road Solon, Ohio USA 44139-4300 440/349-1199FAX: 440/349-1182 www.amresco-inc.com Email: info@amresco-inc.com

CERTIFICATE OF QUALITY

Product: PROTEIN EZ-VISION, 4X

Grade: BIOTECHNOLOGY GRADE

Code: N836

<u>Appearance:</u> Viscous, blue/purple liquid.

Mol. Formula: Not Applicable

Mol. Weight: Not Applicable

Description:

Assay:	<u>Min</u>	<u>Max</u>	<u>Units</u>	Disposition
Abs.@436nm (1:100, Water)	0.2	0.4		PASS
Abs.@560nm (1:100, Water)	0.65	0.96		PASS
Protease		NONE	P/F	PASS

Storage: Store frozen. Keep dry. Warm to room temperature before opening.

n836bio This signature indicates that the above material has met all quality

specifications and has been reviewed by a quality representative.

Signature: You Crocke

Title: QC MANAGER Date: 6/22/2012