



**AMRESKO LLC**

28600 Fountain Parkway  
Solon, Ohio USA 44139  
440/349-1199 FAX: 440/349-1182  
www.amresco-inc.com  
Email: info@amresco-inc.com

ISO 9001 CERTIFIED  
ISO 13485 CERTIFIED

**CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS**

**EZVision ONE, DNA Dye & Buffer, 6X**

**Code: N472**

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable		
Appearance:		Storage:	
Deep red, slightly viscous liquid.		Grade:	BIOTECHNOLOGY GRADE

**Additional Information**

TEST	SPECIFICATION	DISPOSITION
Abs.@340nm (1:100, DI Water)	0.80 - 1.21	PASS
Abs.@521nm (1:500, DI Water)	0.30 - 0.45	PASS
Clarity	PASS	PASS
DNase	NONE	PASS
Electrophoresis (P/F)	PASS	PASS
Protease	NONE	PASS

Spec Set: N472BIO

Analysis may have been rounded to significant digits in specification limits.  
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: \_\_\_\_\_

Title:

Date Printed: 08/12/2016

Page 1 of 1