



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

AMRESCO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Triton X-100

Code: 0694

Chemical Formula:	C34H62O11	Manufacture Date:	(batch specific)
Molecular Weight:	647.00	Expiration/Reassay Date:	(batch specific)
CAS #:	9002-93-1		
Appearance:		Storage:	
Light colored viscous liquid		Grade:	REAGENT GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Cloud Point		PASS
Color (APHA)	<= 100	PASS
Peroxides	NONE	PASS
pH (5%, Water) @25C	6.0 - 8.0	PASS

Spec Set: 0694RG

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed:

04/01/2016

Page 1 of 1