

# SAFETY DATA SHEET

Issuing Date 21-Apr-2009 Revision Date 08-Jan-2016 Revision Number 2

# 1. PRODUCT AND COMPANY IDENTIFICATION

**Product Identifier** 

Product Name Urea

Other means of identification

Product Code(s) 0378

Synonyms No information available

Recommended use of the chemical and restrictions on use

Recommended Use For Further Manufacturing Use Only.
Uses advised against Not for Human or Animal Drug Use

Details of the supplier of the safety data sheet

Company AddressManufacturer AddressVWR International, LLCVWR International, LLCRadnor Corporate Center28600 Fountain Parkway100 Matsonford RoadSolon, Ohio 44139

Radnor, PA 19087-8660

Company Phone Number 1-88-448-4442

E-mail Address info@amresco-inc.com

**Emergency Telephone Number** 

Emergency Telephone Number Chemtrec 1-800-424-9300

# 2. HAZARDS IDENTIFICATION

### Classification

### **OSHA Regulatory Status**

This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.122)

### Label elements

**Emergency Overview** 

Appearance White Physical State Free-flowing Pastille Odor No information available

Hazards not otherwise classified (HNOC)

Not regulated

Other Information

Unknown Acute Toxicity 0 % of the mixture consists of ingredient(s) of unknown toxicity

# 3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS-No	EC No.	Weight %	Trade Secret
Urea	57-13-6	200-315-5	95-100	Not applicable

# 4. FIRST AID MEASURES

**First Aid Measures** 

Eye Contact Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids.

Consult a physician.

**Skin Contact**Wash off immediately with soap and plenty of water removing all contaminated clothes and

shoes.

**Inhalation** Move to fresh air. If not breathing, give artificial respiration.

**Ingestion** Never give anything by mouth to an unconscious person. Rinse mouth.

Most important symptoms and effects, both acute and delayed

Most Important Symptoms/Effects No information available.

Indication of any immediate medical attention and special treatment needed

Notes to Physician Treat symptomatically.

### 5. FIRE-FIGHTING MEASURES

#### **Suitable Extinguishing Media**

Dry chemical, CO<sub>2</sub>, water spray or regular foam.

Unsuitable Extinguishing Media No information available.

### **Specific Hazards Arising from the Chemical**

No information available.

**Hazardous Combustion** 

Carbon oxides. Nitrogen oxides (NOx).

**Products** 

**Explosion Data** 

Sensitivity to Mechanical Impact None. Sensitivity to Static Discharge None.

### **Protective Equipment and Precautions for Firefighters**

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

### 6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

**Personal protection** Avoid contact with skin, eyes and clothing. Ensure adequate ventilation.

**Environmental Precautions** 

**Environmental Precautions** See Section 12 for additional Ecological information.

Methods and material for containment and cleaning up

**Methods for Containment** Prevent further leakage or spillage if safe to do so.

Methods for Cleaning Up Avoid dust formation. Pick up and transfer to properly labeled containers. Ventilate area

and wash spill site after material pickup is complete.

# 7. HANDLING AND STORAGE

No information available

**Precautions for Safe Handling** 

Handling Handle in accordance with good industrial hygiene and safety practice.

Conditions for safe storage, including any incompatibilities

**Storage** Keep containers tightly closed in a dry, cool and well-ventilated place.

Incompatible Products Strong oxidizing agents.

# 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Control parameters

**Exposure Guidelines**This product does not contain any hazardous materials with occupational exposure limits

established by the region specific regulatory bodies.

Chemical Name	ACGIH TLV	OSHA PEL	NIOSH IDLH
Urea	=	-	-
57-13-6			

**Appropriate engineering controls** 

Engineering Measures Showers

Eyewash stations Ventilation systems.

Individual protection measures, such as personal protective equipment

**Eye/Face Protection** Safety glasses with side-shields.

**Skin and Body Protection** Wear protective gloves/clothing.

Respiratory Protection If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved

respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be

provided in accordance with current local regulations.

Hygiene Measures Handle in accordance with good industrial hygiene and safety practice.

# 9. PHYSICAL AND CHEMICAL PROPERTIES

### Information on basic physical and chemical properties

Physical State Free-flowing Pastille

Appearance White Odor

Color No information available Odor Threshold No information available

<u>Property</u> <u>Values</u> <u>Remarks • Method</u>

pH (as 10% solution) 7.2

Melting point/freezing point 132-135 °C

Boiling Point/Range
Flash Point (High in °C)
Evaporation Rate
Flammability (solid, gas)

No information available
No information available
No information available

Flammability Limit in Air

Upper flammability limit: No information available Lower flammability limit: No information available Vapor pressure No information available **Vapor Density** No information available **Specific Gravity** No information available **Water Solubility** Very soluble in water Solubility in other solvents No information available Partition coefficient No information available

### Product Code(s)0378 - Urea

Autoignition TemperatureNo information availableDecomposition TemperatureNo information availableKinematic viscosityNo information availableDynamic viscosityNo information availableExplosive PropertiesNo information availableOxidizing PropertiesNo information available

**Other Information** 

Softening PointNo information availableMolecular WeightNo information availableVOC ContentNo information availableDensityNo information availableBulk DensityNo information available

# 10. STABILITY AND REACTIVITY

#### Reactivity

No data available

#### **Chemical Stability**

Stable under recommended storage conditions.

### **Possibility of Hazardous Reactions**

None under normal processing.

### **Conditions to Avoid**

Extremes of temperature and direct sunlight.

### **Incompatible Materials**

Strong oxidizing agents.

### **Hazardous Decomposition Products**

Carbon oxides. Nitrogen oxides (NOx).

# 11. TOXICOLOGICAL INFORMATION

### Information on likely routes of exposure

#### **Product Information**

**Inhalation** There is no data available for this product.

**Eye Contact** There is no data available for this product.

**Skin Contact** There is no data available for this product.

**Ingestion** No known effect based on information supplied.

Chemical Name	LD50 Oral	LD50 Dermal	LC50 Inhalation
Urea	= 8471 mg/kg (Rat)	-	-
57-13-6	]		

### Information on toxicological effects

Symptoms No information available.

# Delayed and immediate effects as well as chronic effects from short and long-term exposure

SensitizationNo information available.Mutagenic EffectsNo information available.CarcinogenicityNo information available.

Chemical Name	ACGIH	IARC	NTP	OSHA
Urea	-	-	=	-
57-13-6				

Reproductive Toxicity No information available.

STOT - single exposureNo information available.STOT - repeated exposureNo information available.Aspiration hazardNo information available.

Numerical measures of toxicity - Product Information

Unknown Acute Toxicity 0 % of the mixture consists of ingredient(s) of unknown toxicity

The following values are calculated based on chapter 3.1 of the GHS document

# 12. ECOLOGICAL INFORMATION

### **Ecotoxicity**

Chemical Name	Algae/aquatic plants	Toxicity to Fish	Crustacea
Urea	=	16200 - 18300: 96 h Poecilia	10000: 24 h Daphnia magna Straus
57-13-6		reticulata mg/L LC50	mg/L EC50 3910: 48 h Daphnia
		_	magna mg/L EC50 Static

### Persistence and Degradability

No information available.

# Bioaccumulation/Accumulation

No information available.

Chemical Name	Partition coefficient
Urea	-1.59
57-13-6	

Other Adverse Effects No information available

# 13. DISPOSAL CONSIDERATIONS

# Waste treatment methods

Waste Disposal Method Dispose of material in accordance with all federal, state, and local regulations.

**Contaminated Packaging** Do not re-use empty containers.

Chemical Name	RCRA - Halogenated Organic Compounds	RCRA - P Series Wastes	RCRA - F Series Wastes	RCRA - K Series Wastes
Urea 57-13-6	-	-	-	-

Chemical Name	California Hazardous Waste Status
Urea	-
57-13-6	

# 14. TRANSPORT INFORMATION

**DOT** Not regulated

<u>IATA</u> Not regulated

# 15. REGULATORY INFORMATION

**International Inventories** 

TSCA Complies
DSL/NDSL Complies
EINECS/ELINCS Complies

ENCS	Complies
IECSC	Complies
KECL	Complies
PICCS	Complies
AICS	Complies

# Legend:

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

**ENCS** - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

**KECL** - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

AICS - Australian Inventory of Chemical Substances

### U.S. Federal Regulations

#### **SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

Component	SARA 313 - Threshold Values %
Urea	ē
57-13-6 ( 95-100 )	

# SARA 311/312 Hazard Categories

Acute Health HazardNoChronic Health HazardNoFire HazardNoSudden Release of Pressure HazardNoReactive HazardNo

#### **Clean Water Act**

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

Chemical Name	CWA - Reportable Quantities	CWA - Toxic Pollutants	CWA - Priority Pollutants	CWA - Hazardous Substances
Urea 57-13-6	-	-	-	-

# CERCLA

This material, as supplied, does not contain any substances regulated as hazardous substances under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40 CFR 355). There may be specific reporting requirements at the local, regional, or state level pertaining to releases of this material

Chemical Name	Hazardous Substances RQs	Extremely Hazardous Substances RQs	RQ
Urea 57-13-6	-	-	•

# U.S. State Regulations

# **California Proposition 65**

This product does not contain any Proposition 65 chemicals.

Chemical Name	California Prop. 65	
Urea - 57-13-6	-	

### U.S. State Right-to-Know Regulations

Chemical Name	New Jersey	Massachusetts	Pennsylvania
Urea	-	-	-
57-13-6			

### U.S. EPA Label Information

# EPA Pesticide Registration Number Not regulated

# **16. OTHER INFORMATION**

Issuing Date21-Apr-2009Revision Date08-Jan-2016

**Revision Note** 

No information available

Disclaimer

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. VWR International and its Affiliates shall not be held liable for any damage resulting from handling.

**End of Safety Data Sheet**