



AMRESCO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Urea

Code: 0378

Chemical Formula:	CH4N2O	Manufacture Date:	(batch specific)
Molecular Weight:	60.10	Expiration/Reassay Date:	(batch specific)
CAS #:	57-13-6		
Appearance:		Storage:	Room Temperature
Fine, white, free-flowing pastilles.		Grade:	USP GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Chloride	<= 0.007 %	PASS
Heavy Metals (as Pb)	<= 0.002 %	PASS
Identification	PASS	PASS
Insolubles (Alcohol)	<= 0.04 %	PASS
Melting Range	132 - 135 Degrees C	PASS
Purity	99.0 - 100.5 %	PASS
Residue after Ignition	<= 0.1 %	PASS
Sulfate	<= 0.010 %	PASS

Spec Set: 0378USP

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Residual Solvents, if applicable will be reported for actual batch number only.

Internal ID #: 269

Signature: _____
Title:

Date Printed: 03/09/2016
Page 1 of 1