



AMRESKO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Polyethylene Glycol 8000

Code: 0159

| | | | |
|--------------------|---|---------------------------------|---------------------|
| Chemical Formula: | HOCH ₂ CH ₂ (OCH ₂ CH ₂) _n OH | Manufacture Date: | (batch specific) |
| Molecular Weight: | | Expiration/Reassay Date: | (batch specific) |
| CAS #: | 25322-68-3 | Storage: | |
| Appearance: | | Grade: | BIOTECHNOLOGY GRADE |
| White waxy powder. | | Additional Product Description: | |
| | | Molecular Weight range is | 7000-9000. |

Additional Information

| TEST | SPECIFICATION | DISPOSITION |
|----------------------|-------------------|-------------|
| DNase | NONE | PASS |
| pH (5%, Water) @25C | 4.5 - 7.5 | PASS |
| RNase | NONE | PASS |
| Solidification Point | 55 - 64 Degrees C | PASS |

Spec Set: 0159BIO

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____
Title:

Date Printed: 02/29/2016
Page 1 of 1