



**AMRESKO LLC**

28600 Fountain Parkway  
Solon, Ohio USA 44139  
440/349-1199 FAX: 440/349-1182  
www.amresco-inc.com  
Email: info@amresco-inc.com

ISO 9001 CERTIFIED  
ISO 13485 CERTIFIED

**CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS**

**TEMED**

**Code: 0761**

Chemical Formula:	C6H16N2	Manufacture Date:	(batch specific)
Molecular Weight:	116.21	Expiration/Reassay Date:	(batch specific)
CAS #:	110-18-9		
Appearance:		Storage:	
Clear liquid		Grade:	ULTRA PURE GRADE

**Additional Information**

TEST	SPECIFICATION	DISPOSITION
Boiling Range	119 - 121 Degrees C	PASS
Purity (Anhydrous Basis)	>= 99.0 %	PASS
Refractive Index	1.417 - 1.419	PASS

Spec Set: 0761UP

Analysis may have been rounded to significant digits in specification limits.  
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: \_\_\_\_\_  
Title:

Date Printed: 06/07/2016  
Page 1 of 1