



AMRESCO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

BCIP/NBT AMP Buffer Solution

Code: E116

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable	Storage:	
Appearance:		Grade:	BIOTECHNOLOGY GRADE
Clear yellow liquid		Additional Product Description:	
			0.42 g/L NBT, and 0.21 g/L BCIP in 0.1M AMP Buffer Solution.

Additional Information

TEST	SPECIFICATION	DISPOSITION
Abs.@234nm (3%, Water)	1.07 - 1.37	PASS
Abs.@260nm (3%, Water)	0.86 - 1.09	PASS
Conductivity	2500 - 5500 umhos	PASS
Moisture (by Karl Fisher)		PASS
pH @25C	9.7 - 9.9	PASS

Spec Set: E116BIO

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____
Title:

Date Printed: 06/09/2016
Page 1 of 1