



**AMRESKO LLC**

28600 Fountain Parkway  
Solon, Ohio USA 44139  
440/349-1199 FAX: 440/349-1182  
www.amresco-inc.com  
Email: info@amresco-inc.com

ISO 9001 CERTIFIED  
ISO 13485 CERTIFIED

---

**CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS**

**Acryl/Bis Solution (30%) 37.5:1**

**Code: E347**

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable	Storage:	
Appearance:		Grade:	ULTRA PURE GRADE
Clear and colorless solution		Additional Product Description:	
			A 30% (W/V) solution of acrylamide (29.22%) and bis-acrylamide (0.78%) in specially prepared deionized water. Final ratio is 37.5:1.

---

**Additional Information**

TEST	SPECIFICATION	DISPOSITION
Conductivity	<= 10 umhos	PASS
pH @25C		PASS

Spec Set: E347UP

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: \_\_\_\_\_

Title:

Date Printed:

06/10/2016

Page 1 of 1