



**AMRESKO LLC**

28600 Fountain Parkway  
Solon, Ohio USA 44139  
440/349-1199 FAX: 440/349-1182  
www.amresco-inc.com  
Email: info@amresco-inc.com

ISO 9001 CERTIFIED  
ISO 13485 CERTIFIED

---

**CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS**

**SDS 20% SOLUTION**

**Code: 0837**

Chemical Formula:	C <sub>12</sub> H <sub>25</sub> NaO <sub>4</sub> S	Manufacture Date:	(batch specific)
Molecular Weight:	288.40	Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable		
Appearance:		Storage:	
Clear, particulate-free, colorless liquid		Grade:	BIOTECHNOLOGY GRADE
		Additional Product Description:	
			A solution that is 20% Sodium Dodecyl Sulfate. MW and MF noted are for the active ingredient.

---

**Additional Information**

TEST	SPECIFICATION	DISPOSITION
Conductivity	22800 - 27600 umhos	PASS
pH @25C	5.0 - 8.0	PASS

Spec Set: 0837BIO

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: \_\_\_\_\_

Title:

Date Printed:

06/08/2016

Page 1 of 1