



ISO 9001 CERTIFIED  
ISO 13485 CERTIFIED

## AMRESCO LLC

28600 Fountain Parkway  
Solon, Ohio USA 44139  
440/349-1199 FAX: 440/349-1182  
www.amresco-inc.com  
Email: info@amresco-inc.com

### CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

#### Yeast Nitrogen

**Code:** J630

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable		
Appearance:		Storage:	
Light yellow to beige powder		Grade:	BACTERIOLOGICAL

#### Additional Information

TEST	SPECIFICATION	DISPOSITION
Loss on Drying	$\leq 4.0 \%$	PASS
pH (1:10, 1.7g/100ml Water) 25C	4.3 - 4.7	PASS
Solubility (1.7g/100ml Water)	PASS	PASS

Spec Set: J630BACT

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: \_\_\_\_\_

Title:

Date Printed:

03/07/2016

Page 1 of 1