



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

AMRESO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Ammonium Persulfate

Code: 0486

Chemical Formula:	(NH ₄) ₂ S ₂ O ₈	Manufacture Date:	(batch specific)
Molecular Weight:	228.19	Expiration/Reassay Date:	(batch specific)
CAS #:	7727-54-0		
Appearance:		Storage:	
White crystalline powder		Grade:	ACS GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Chloride & Chlorate	<= 0.001 %	PASS
Heavy Metals (as Pb)	<= 0.005 %	PASS
Insolubles	<= 0.005 %	PASS
Iron	<= 0.001 %	PASS
Manganese	<= 0.00005 %	PASS
Purity	>= 98.0 %	PASS
Residue on Ignition	<= 0.05 %	PASS
Titration Free Acid	<= 0.04 meq/g	PASS

Spec Set: 0486ACS

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed:

05/04/2016

Page 1 of 1