



ISO 9001 CERTIFIED  
ISO 13485 CERTIFIED

## AMRESKO LLC

28600 Fountain Parkway  
Solon, Ohio USA 44139  
440/349-1199 FAX: 440/349-1182  
www.amresco-inc.com  
Email: info@amresco-inc.com

### CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

#### Methanol Anhydrous

**Code:** 0323

Chemical Formula:	CH <sub>3</sub> OH	Manufacture Date:	(batch specific)
Molecular Weight:	32.04	Expiration/Reassay Date:	(batch specific)
CAS #:	67-56-1		
Appearance:		Storage:	
Clear, colorless liquid		Grade:	ULTRA PURE GRADE

#### Additional Information

TEST	SPECIFICATION	DISPOSITION
Moisture (KF)	$\leq 0.1 \%$	PASS
Purity	$\geq 99.8 \%$	PASS
Specific Gravity	0.783 - 0.791	PASS

Spec Set: 0323UP

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: \_\_\_\_\_

Title:

Date Printed:

03/07/2016

Page 1 of 1