



AMRESKO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Puromycin Dihydrochloride

Code: J593

Chemical Formula:	C22H29N7O5.2HCl	Manufacture Date:	(batch specific)
Molecular Weight:	544.40	Expiration/Reassay Date:	(batch specific)
CAS #:	58-58-2		
Appearance:		Storage:	
White to yellow-white solid		Grade:	ULTRA PURE GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Melting Point	168 - 178 Degrees C	PASS
Purity	>= 98 %	PASS
Solubility (5%, Water)	PASS	PASS

Spec Set: J593UP

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____
Title:

Date Printed: 06/23/2016
Page 1 of 1