



AMRESKO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Sodium Phosphate, Dibasic, Anhydrous

Code: 0404

Chemical Formula:	Na2HPO4	Manufacture Date:	(batch specific)
Molecular Weight:	141.96	Expiration/Reassay Date:	(batch specific)
CAS #:	7558-79-4		
Appearance:		Storage:	Room Temperature
White powder		Grade:	ACS GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Chloride	<= 0.002 %	PASS
Heavy Metals (as Pb)	<= 0.001 %	PASS
Insolubles	<= 0.01 %	PASS
Iron	<= 0.002 %	PASS
Loss on Drying @105C	<= 0.2 %	PASS
pH (5%, Water) @25C	8.7 - 9.3	PASS
Purity	>= 99.0 %	PASS
Sulfate	<= 0.005 %	PASS

Spec Set: 0404ACS

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed: 05/04/2016

Page 1 of 1