



AMRESKO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

SDS 20% SOLUTION

Code: M112

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable	Storage:	
Appearance:		Grade:	PROTEOMICS GRADE
Clear, particulate-free, colorless liquid		Additional Product Description:	
		A solution that is 20% Sodium Dodecyl Sulfate.	

Additional Information

TEST	SPECIFICATION	DISPOSITION
DNase	NONE	PASS
pH @25C	5.0 - 8.0	PASS
Protease	NONE	PASS
RNase	NONE	PASS

Spec Set: M112PTM

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____
Title:

Date Printed: 08/10/2016
Page 1 of 1