



AMRESKO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

TBE Dry Powder Buffer

Code: 0478

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable		
Appearance:		Storage:	
White powder		Grade:	ULTRA PURE GRADE
		Additional Product Description:	
			A single strength (1X) solution contains 0.089 M Tris base, 0.089 M borate & 0.002 M EDTA.

Additional Information

TEST	SPECIFICATION	DISPOSITION
DNase	NONE	PASS
pH (1X Concentration) @25C	8.1 - 8.5	PASS
Protease	NONE	PASS
RNase	NONE	PASS
Solubility (1X Concentration)	PASS	PASS

Spec Set: 0478UP

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed: 03/08/2016

Page 1 of 1