



ISO 9001 CERTIFIED  
ISO 13485 CERTIFIED

## AMRESKO LLC

28600 Fountain Parkway  
Solon, Ohio USA 44139  
440/349-1199 FAX: 440/349-1182  
www.amresco-inc.com  
Email: info@amresco-inc.com

### CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

#### Potassium Acetate

**Code:** 0698

Chemical Formula:	C <sub>2</sub> H <sub>3</sub> O <sub>2</sub> K	Manufacture Date:	(batch specific)
Molecular Weight:	98.14	Expiration/Reassay Date:	(batch specific)
CAS #:	127-08-2		
Appearance:		Storage:	
White crystalline powder		Grade:	ACS GRADE

#### Additional Information

TEST	SPECIFICATION	DISPOSITION
Calcium	<= 0.005 %	PASS
Chloride	<= 0.003 %	PASS
Heavy Metals	<= 0.0005 %	PASS
Insolubles	<= 0.005 %	PASS
Iron	<= 0.0005 %	PASS
Magnesium	<= 0.002 %	PASS
pH (5%, Water) @25C	6.5 - 9.0	PASS
Phosphate	<= 0.001 %	PASS
Purity	>= 99.0 %	PASS
Sodium	<= 0.03 %	PASS
Sulfate	<= 0.002 %	PASS

Spec Set: 0698ACS

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: \_\_\_\_\_

Title:

Date Printed:

03/08/2016

Page 1 of 1