



AMRESKO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Ponceau S Stain

Code: K793

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable	Storage:	
Appearance:		Grade:	PROTEOMICS GRADE
Red, haze-free liquid.			

Additional Information

TEST	SPECIFICATION	DISPOSITION
Abs.@ Lambda Max (1:50, Water)	0.9000 - 1.4000	PASS
Lambda Max (1:50, Water)	517 - 523 nm	PASS
pH (1:50, Water) @25C		PASS

Spec Set: K793PTM

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____
Title:

Date Printed: 07/29/2016
Page 1 of 1