



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

AMRESCO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

RiboZol, RNA Extraction Reagent

Code: N580

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable		
Appearance:		Storage:	Room Temperature
Pink to dark pink, clear liquid.		Grade:	BIOTECHNOLOGY GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Abs.@ Lambda Max	1.50 - 2.70	PASS
DNase	NONE	PASS
Electrophoresis (P/F)	PASS	PASS
Lambda Max	530 - 534 nm	PASS
Percent OD260/280	>= 95 %	PASS
Percent Yield	>= 80 %	PASS
pH @25C	3.8 - 4.2	PASS
RNase	NONE	PASS
Specific Gravity @20C	1.0824 - 1.1124	PASS
Water Content (KF)	36 - 42 %	PASS

Spec Set: N580BIO

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed:

08/12/2016

Page 1 of 1