



DISPOSABLE HOMOGENIZER BioMasher®

Made in Japan


BioMasher-I





3 The pestle is necessary.




BioMasher-II



BioMasher-III



BioMasher-IV



	BioMasher-I	BioMasher-II	BioMasher-III	BioMasher-IV
Description	For many samples	For increase of the collection amount	For increase of the collection amount For many samples	For hard samples For many samples
Collection amount	Small	Large	Large	Large
Heat	Heatless	Controllable	Controllable	Heat
Device	Centrifuge	Unnecessary	Centrifuge	Cell crusher
Soft samples	◎	◎	◎	○
Hard samples	△	◎	◎	◎
Yeast	×	×	×	◎
Cultured cell	△※	×	△	○

※Please inquire for the usage.

		Cat.No.	Sets/ Pack
BioMasher I for Soft samples			
1.5ml Microtube, with O-ring, PP filter		NIP-30-1.5-O	30
		NIP-120-1.5-O	120
2.0ml Microtube, with O-ring, PP filter		NIP-30-2.0-O	30
		NIP-120-2.0-O	120
BioMasher I for Hard samples			
1.5ml Microtube, without O-ring, PP filter		NIP-30-1.5	30
		NIP-120-1.5	120
2.0ml Microtube, without O-ring, PP filter		NIP-30-2.0	30
		NIP-120-2.0	120
BioMasher I with small pore size filter for Soft samples			
1.5ml Microtube, with O-ring, PE filter		NIP-30-1.5-OS	30
		NIP-120-1.5-OS	120
2.0ml Microtube, with O-ring, PE filter		NIP-30-2.0-OS	30
		NIP-120-2.0-OS	120
BioMasher I with small pore size filter for Hard samples			
1.5ml Microtube, without O-ring, PE filter		NIP-30-1.5-NS	30
		NIP-120-1.5-NS	120
2.0ml Microtube, without O-ring, PE filter		NIP-30-2.0-NS	30
		NIP-120-2.0-NS	120
BioMasher II	Non-sterile	320102	100
	EOG-sterilized	320103	100
BioMasher III	Non-sterile	320302	50
BioMasher IV	Non-sterile	320402	240
	EOG-sterilized	320403	240
BioMasher V	Non-sterile	891390	20
		891391	100
	EOG-sterilized	891392	20
		891393	100
PowerMasher	Pestle motor for BioMasher II, III & V	320105	1pcs
PowerMasher II	Pestle motor for BioMasher II, III & V	891300	1pcs

Frontiers in Life Science

Your Local Dealer

Distributed by

 **funakoshi**
FRONTIERS IN LIFE SCIENCE

URL : <http://www.funakoshi.co.jp/export/index.php>

e-mail : export@funakoshi.co.jp

Tel. : +81-3-5684-6296, Fax : +81-3-5684-6297