

**AMRESCO LLC**

28600 Fountain Parkway
 Solon, Ohio USA 44139
 440/349-1199 FAX: 440/349-1182
 www.amresco-inc.com
 Email: info@amresco-inc.com

ISO 9001 CERTIFIED
 ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS**Ammonium Sulfate****Code: 0191**

Chemical Formula:	H8N2O4S	Manufacture Date:	(batch specific)
Molecular Weight:	132.14	Expiration/Reassay Date:	(batch specific)
CAS #:	7783-20-2		
Appearance:		Storage:	Room Temperature
White crystalline powder		Grade:	ULTRA PURE GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Arsenic	<= 0.0001 %	PASS
Chloride	<= 0.0005 %	PASS
DNase	NONE	PASS
Heavy Metals (as Pb)	<= 0.0002 %	PASS
Insolubles	<= 0.005 %	PASS
Iron	<= 0.0005 %	PASS
Magnesium Test	PASS	PASS
pH (5%, Water) @25C	5.0 - 6.0	PASS
Phosphate	<= 0.0005 %	PASS
Protease	NONE	PASS
Purity (Dry Basis)	>= 99.5 %	PASS
Residue on Ignition	<= 0.005 %	PASS
RNase	NONE	PASS
Selenium	<= 0.0020 %	PASS

Spec Set: 0191UP

Analysis may have been rounded to significant digits in specification limits.
 Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed: 02/29/2016

Page 1 of 1