



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

AMRESCO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

CHAPS

Code: 0465

Chemical Formula:	C32H58N2SO7	Manufacture Date:	(batch specific)
Molecular Weight:	614.89	Expiration/Reassay Date:	(batch specific)
CAS #:	75621-03-3		
Appearance:		Storage:	Room Temperature
White crystalline powder		Grade:	ULTRA PURE GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Conductivity (10%, Water)	<= 50 umhos	PASS
DNase	NONE	PASS
Identification	PASS	PASS
pH (10%, Water) @25C	5.0 - 7.0	PASS
Protease	NONE	PASS
Purity (via TLC)	ONE SPOT	PASS
Residue on Ignition	<= 0.1 %	PASS
RNase	NONE	PASS
Solubility (20%, Water)	PASS	PASS

Spec Set: 0465UP

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed:

03/08/2016

Page 1 of 1