



AMRESCO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Hygromycin B

Code: J607

Chemical Formula:	C20H37N3O13	Manufacture Date:	(batch specific)
Molecular Weight:	527.50	Expiration/Reassay Date:	(batch specific)
CAS #:	31282-04-9		
Appearance:		Storage:	Cold
White to pale-brown powder		Grade:	ULTRA PURE GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Activity	>= 1000 U/mg	PASS
Purity (HPLC)	>= 80 %	PASS

Spec Set: J607UP

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____
Title:

Date Printed: 06/23/2016
Page 1 of 1