



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

AMRESCO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

POTASSIUM PHOSPHATE DIBASIC ANHYDROUS

Code: 0705

Chemical Formula:	K ₂ HPO ₄	Manufacture Date:	(batch specific)
Molecular Weight:	174.18	Expiration/Reassay Date:	(batch specific)
CAS #:	7758-11-4		
Appearance:		Storage:	Room Temperature
Fine white crystalline powder		Grade:	ACS GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Chloride	<= 0.003 %	PASS
Heavy Metals (as Pb)	<= 0.0005 %	PASS
Insolubles	<= 0.01 %	PASS
Iron	<= 0.001 %	PASS
Loss on Drying	<= 1.0 %	PASS
Nitrogen Compounds	<= 0.001 %	PASS
pH (5%, Water) @25C	8.5 - 9.6	PASS
Purity	>= 98.0 %	PASS
Sodium	<= 0.05 %	PASS
Sulfate	<= 0.005 %	PASS

Spec Set: 0705ACS

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed:

03/08/2016

Page 1 of 1