



AMRESCO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

NEXT GEL, 7.5%

Code: M255

Chemical Formula:	Not Applicable	Manufacture Date:	(batch specific)
Molecular Weight:		Expiration/Reassay Date:	(batch specific)
CAS #:	Not Applicable		
Appearance:		Storage:	
Clear, colorless liquid.		Grade:	PROTEOMICS GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Conductivity (1:20, Water)	430 - 550 uS	PASS
Electrophoresis (P/F)	PASS	PASS
pH @25C	6.8 - 7.2	PASS
Polymerization Time	<= 30 minutes	PASS
Protease	NONE	PASS

Spec Set: M255PTM

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____
Title:

Date Printed: 08/11/2016
Page 1 of 1