



ISO 9001 CERTIFIED
 ISO 13485 CERTIFIED

AMRESCO LLC
 28600 Fountain Parkway
 Solon, Ohio USA 44139
 440/349-1199 FAX: 440/349-1182
 www.amresco-inc.com
 Email: info@amresco-inc.com

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

POTASSIUM PHOSPHATE, MONOBASIC, ANHYDROUS

Code: 0781

Chemical Formula:	KH ₂ PO ₄	Manufacture Date:	(batch specific)
Molecular Weight:	136.09	Expiration/Reassay Date:	(batch specific)
CAS #:	7778-77-0		
Appearance:		Storage:	Room Temperature
White solid		Grade:	ACS GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Chloride	≤ 0.001 %	PASS
Heavy Metals	≤ 0.001 %	PASS
Insolubles	≤ 0.01 %	PASS
Iron	≤ 0.002 %	PASS
Loss on Drying	≤ 0.2 %	PASS
Nitrogen Compounds	≤ 0.001 %	PASS
pH (5%, Water) @25C	4.1 - 4.5	PASS
Purity	≥ 99.0 %	PASS
Sodium	≤ 0.005 %	PASS
Sulfate	≤ 0.003 %	PASS

Spec Set: 0781ACS

Analysis may have been rounded to significant digits in specification limits.
 Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____
 Title:

Date Printed: 06/07/2016
 Page 1 of 1