

## **SAFETY DATA SHEET**

Issuing Date 07-Jun-2010 Revision Date 08-Jan-2016 Revision Number 2

## 1. PRODUCT AND COMPANY IDENTIFICATION

**Product Identifier** 

Product Name Ethidium Bromide 0.625 mg/mL

Other means of identification

Product Code(s) E406

Synonyms No information available

Recommended use of the chemical and restrictions on use

Recommended Use For Further Manufacturing Use Only.
Uses advised against Not for Human or Animal Drug Use

Details of the supplier of the safety data sheet

Company AddressManufacturer AddressVWR International, LLCVWR International, LLCRadnor Corporate Center28600 Fountain Parkway100 Matsonford RoadSolon, Ohio 44139

Radnor, PA 19087-8660

Company Phone Number1-800-448-4442E-mail Addressinfo@amresco-inc.com

**Emergency Telephone Number** 

Emergency Telephone Number Chemtrec 1-800-424-9300

## 2. HAZARDS IDENTIFICATION

#### Classification

**OSHA Regulatory Status** 

This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.122)

## Label elements

**Emergency Overview** 

Appearance Dark red Physical State Liquid Odor No information available

Hazards not otherwise classified (HNOC)

Not regulated

Other Information

## 3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS-No	EC No.	Weight %	Trade Secret
3,8-Diamino-1-ethyl-6-ph	1239-45-8	214-984-6	<0.1	Not applicable
enylphenantridinium				
bromide				

#### 4. FIRST AID MEASURES

**First Aid Measures** 

Eye Contact Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids.

Consult a physician.

**Skin Contact**Wash off immediately with soap and plenty of water removing all contaminated clothes and

shoes

**Inhalation** Move to fresh air. If not breathing, give artificial respiration.

**Ingestion** Never give anything by mouth to an unconscious person. Rinse mouth.

Most important symptoms and effects, both acute and delayed

Most Important Symptoms/Effects No information available.

Indication of any immediate medical attention and special treatment needed

Notes to Physician Treat symptomatically.

#### 5. FIRE-FIGHTING MEASURES

#### **Suitable Extinguishing Media**

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable Extinguishing Media No information available.

## **Specific Hazards Arising from the Chemical**

No information available.

**Explosion Data** 

Sensitivity to Mechanical Impact None. Sensitivity to Static Discharge None.

## **Protective Equipment and Precautions for Firefighters**

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

#### 6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

**Personal protection** Ensure adequate ventilation, especially in confined areas.

**Environmental Precautions** 

**Environmental Precautions** See Section 12 for additional Ecological information.

Methods and material for containment and cleaning up

**Methods for Containment** Prevent further leakage or spillage if safe to do so.

Methods for Cleaning Up Absorb with sand or vermiculite. Pick up and transfer to properly labeled containers.

Ventilate area and wash spill site after material pickup is complete.

## 7. HANDLING AND STORAGE

#### **Precautions for Safe Handling**

## Product Code(s)E406 - Ethidium Bromide 0.625 mg/mL

Handling Handle in accordance with good industrial hygiene and safety practice.

Conditions for safe storage, including any incompatibilities

**Storage** Keep containers tightly closed in a dry, cool and well-ventilated place.

Incompatible Products Strong oxidizing agents.

#### 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Control parameters

Exposure Guidelines This product does not contain any hazardous materials with occupational exposure limits

established by the region specific regulatory bodies.

Chemical Name	ACGIH TLV	OSHA PEL	NIOSH IDLH
3,8-Diamino-1-ethyl-6-phenylphena	-	=	-
ntridinium bromide			
1239-45-8			

**Appropriate engineering controls** 

Engineering Measures Showers

Eyewash stations Ventilation systems.

Individual protection measures, such as personal protective equipment

**Eye/Face Protection** Avoid contact with eyes. Safety glasses with side-shields.

**Skin and Body Protection** Wear protective gloves/clothing.

Respiratory Protection If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved

respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be

Remarks • Method

provided in accordance with current local regulations.

**Hygiene Measures** Handle in accordance with good industrial hygiene and safety practice.

## 9. PHYSICAL AND CHEMICAL PROPERTIES

#### Information on basic physical and chemical properties

Physical State Liquid

AppearanceDark redOdorNo information availableColorNo information availableOdor ThresholdNo information available

Property Values

pH No information available

Melting point/freezing point

Boiling Point/Range No information available Flash Point (High in °C) No information available Evaporation Rate No information available Flammability (solid, gas) No information available

Flammability Limit in Air Upper flammability limit:

Upper flammability limit:
Lower flammability limit:
Vapor pressure
Vapor Density

No information available
No information available
No information available

Specific Gravity 0.9995

Water Solubility
Solubility in other solvents
No information available
No information available

## Product Code(s)E406 - Ethidium Bromide 0.625 mg/mL

Partition coefficient
Autoignition Temperature
Decomposition Temperature
Kinematic viscosity
Dynamic viscosity
Explosive Properties
Oxidizing Properties
No information available

Other Information

Softening PointNo information availableMolecular WeightNo information availableVOC ContentNo information availableDensityNo information availableBulk DensityNo information available

## 10. STABILITY AND REACTIVITY

## Reactivity

No data available

#### **Chemical Stability**

Stable under recommended storage conditions.

## **Possibility of Hazardous Reactions**

None under normal processing.

## **Conditions to Avoid**

Extremes of temperature and direct sunlight.

## **Incompatible Materials**

Strong oxidizing agents.

#### **Hazardous Decomposition Products**

Carbon oxides. Nitrogen oxides (NOx). Hydrogen bromide gas.

## 11. TOXICOLOGICAL INFORMATION

#### Information on likely routes of exposure

#### **Product Information**

**Inhalation** There is no data available for this product.

**Eye Contact** There is no data available for this product.

**Skin Contact** There is no data available for this product.

**Ingestion** There is no data available for this product.

Chemical Name	LD50 Oral	LD50 Dermal	LC50 Inhalation
3,8-Diamino-1-ethyl-6-phenylphena	-	-	-
ntridinium bromide			
1239-45-8			

#### Information on toxicological effects

**Symptoms** No information available.

## Delayed and immediate effects as well as chronic effects from short and long-term exposure

SensitizationNo information available.Mutagenic EffectsNo information available.CarcinogenicityNo information available.

# Product Code(s)E406 - Ethidium Bromide 0.625 mg/mL

Chemical Name	ACGIH	IARC	NTP	OSHA
3,8-Diamino-1-ethyl-6-pheny	-	-	-	-
lphenantridinium bromide				
1239-45-8				

Reproductive Toxicity
STOT - single exposure
STOT - repeated exposure
Aspiration hazard
No information available.
No information available.
No information available.

Numerical measures of toxicity - Product Information

## 12. ECOLOGICAL INFORMATION

#### **Ecotoxicity**

100% of the mixture consists of components(s) of unknown hazards to the aquatic environment

Chemical Name	Algae/aquatic plants	Toxicity to Fish	Crustacea
3,8-Diamino-1-ethyl-6-phenylphena	-	-	-
ntridinium bromide			
1239-45-8			

#### Persistence and Degradability

No information available.

#### **Bioaccumulation/Accumulation**

No information available.

Chemical Name	Partition coefficient
3,8-Diamino-1-ethyl-6-phenylphenantridinium bromide 1239-45-8	-

Other Adverse Effects No information available

## 13. DISPOSAL CONSIDERATIONS

#### Waste treatment methods

Waste Disposal Method Dispose of material in accordance with all federal, state, and local regulations.

**Contaminated Packaging** Do not re-use empty containers.

Chemical Name	RCRA - Halogenated Organic Compounds	RCRA - P Series Wastes	RCRA - F Series Wastes	RCRA - K Series Wastes
3,8-Diamino-1-ethyl-6-pheny	-	-	=	-
Iphenantridinium bromide				
1239-45-8				

Chemical Name	California Hazardous Waste Status
3,8-Diamino-1-ethyl-6-phenylphenantridinium bromide	-
1239-45-8	

## 14. TRANSPORT INFORMATION

DOTNot regulatedIATANot regulated

## **15. REGULATORY INFORMATION**

## Product Code(s)E406 - Ethidium Bromide 0.625 mg/mL

#### **International Inventories**

**TSCA** Does not comply **DSL/NDSL** Does not comply Complies **EINECS/ELINCS** Does not Comply **ENCS** Complies **IECSC** Complies **KECL** Complies **PICCS** Complies **AICS** 

#### Legend:

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

**ENCS** - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

**KECL** - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

AICS - Australian Inventory of Chemical Substances

#### U.S. Federal Regulations

#### **SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

Component	SARA 313 - Threshold Values %
3,8-Diamino-1-ethyl-6-phenylphenantridinium bromide	-
1239-45-8 (<0.1)	

#### SARA 311/312 Hazard Categories

Acute Health Hazard	No
Chronic Health Hazard	No
Fire Hazard	No
Sudden Release of Pressure Hazard	No
Reactive Hazard	No

## **Clean Water Act**

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

Chemical Name	CWA - Reportable Quantities	CWA - Toxic Pollutants	CWA - Priority Pollutants	CWA - Hazardous Substances
3,8-Diamino-1-ethyl-6-pheny lphenantridinium bromide	-	-	-	-
1239-45-8				

#### **CERCLA**

This material, as supplied, does not contain any substances regulated as hazardous substances under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40 CFR 355). There may be specific reporting requirements at the local, regional, or state level pertaining to releases of this material

Chemical Name	Hazardous Substances RQs	Extremely Hazardous Substances RQs	RQ
3,8-Diamino-1-ethyl-6-phenylphena	=	-	-
ntridinium bromide			
1239-45-8			

#### **U.S. State Regulations**

#### **California Proposition 65**

This product does not contain any Proposition 65 chemicals.

Chemical Name	California Prop. 65	

_		
I	3,8-Diamino-1-ethyl-6-phenylphenantridinium bromide - 1239-45-8	_
- 1	3,0-Diamino-1-emyl-o-phenylphenanthum brothide - 1233-43-0	

#### **U.S. State Right-to-Know Regulations**

Chemical Name	New Jersey	Massachusetts	Pennsylvania
3,8-Diamino-1-ethyl-6-phenylphena	-	-	-
ntridinium bromide			
1239-45-8			

## U.S. EPA Label Information

EPA Pesticide Registration Number Not regulated

## **16. OTHER INFORMATION**

Issuing Date07-Jun-2010Revision Date08-Jan-2016

Revision Note

No information available

**Disclaimer** 

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. VWR International and its Affiliates shall not be held liable for any damage resulting from handling.

**End of Safety Data Sheet**