



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

AMRESCO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Trichloroacetic Acid

Code: 0133

Chemical Formula:	CCl ₃ COOH	Manufacture Date:	(batch specific)
Molecular Weight:	163.39	Expiration/Reassay Date:	(batch specific)
CAS #:	76-03-9		
Appearance:		Storage:	Cold
White deliquescent crystals		Grade:	ACS GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Chloride	<= 0.002 %	PASS
Clarity of Solution	PASS	PASS
Heavy Metals	<= 0.002 %	PASS
Insoluble Matter	<= 0.01 %	PASS
Iron	<= 0.001 %	PASS
Nitrate	<= 0.002 %	PASS
Phosphate	<= 5 ppm	PASS
Purity	>= 99.0 %	PASS
Residue after Ignition	<= 0.03 %	PASS
Substances Darkened by H ₂ SO ₄	PASS	PASS
Sulfate	<= 0.02 %	PASS

Spec Set: 0133ACS

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed:

02/26/2016

Page 1 of 1