



**AMRESCO LLC**

28600 Fountain Parkway  
Solon, Ohio USA 44139  
440/349-1199 FAX: 440/349-1182  
www.amresco-inc.com  
Email: info@amresco-inc.com

ISO 9001 CERTIFIED  
ISO 13485 CERTIFIED

**CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS**

**Sodium Azide**

**Code: 0639**

Chemical Formula:	Na N <sub>3</sub>	Manufacture Date:	(batch specific)
Molecular Weight:	65.01	Expiration/Reassay Date:	(batch specific)
CAS #:	26628-22-8		
Appearance:		Storage:	
White crystalline powder		Grade:	HIGH PURITY GRADE

**Additional Information**

TEST	SPECIFICATION	DISPOSITION
DNase	NONE	PASS
Heavy Metals	<= 0.002 %	PASS
Iron	<= 0.001 %	PASS
Loss on Drying	<= 0.5 %	PASS
Melting Range w/Decomposition	>= 320 Degrees C	PASS
Protease	NONE	PASS
Purity	>= 99.0 %	PASS
RNase	NONE	PASS
Solubility (20%, Water)	PASS	PASS

Spec Set: 0639HPG

Analysis may have been rounded to significant digits in specification limits.  
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: \_\_\_\_\_  
Title:

Date Printed: 05/12/2016  
Page 1 of 1