



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

AMRESCO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Boric Acid

Code: 0588

Chemical Formula:	H ₃ BO ₃	Manufacture Date:	(batch specific)
Molecular Weight:	61.83	Expiration/Reassay Date:	(batch specific)
CAS #:	10043-35-3		
Appearance:		Storage:	
White granular powder		Grade:	ACS GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Calcium	<= 0.005 %	PASS
Chloride	<= 0.001 %	PASS
Heavy Metals (as Pb)	<= 0.001 %	PASS
Insolubles (Methanol)	<= 0.005 %	PASS
Iron	<= 0.001 %	PASS
Nonvolatiles (Methanol)	<= 0.05 %	PASS
Phosphate	<= 0.001 %	PASS
Purity	>= 99.5 %	PASS
Sulfate	<= 0.01 %	PASS

Spec Set: 0588ACS

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____

Title:

Date Printed:

05/04/2016

Page 1 of 1