



**AMRESKO LLC**

28600 Fountain Parkway  
Solon, Ohio USA 44139  
440/349-1199 FAX: 440/349-1182  
www.amresco-inc.com  
Email: info@amresco-inc.com

ISO 9001 CERTIFIED  
ISO 13485 CERTIFIED

---

**CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS**

**5/175 mg pNPP Substrate Tabs**

**Code: 0405**

Chemical Formula:	C6H4NO6PNa2.6H2O	Manufacture Date:	(batch specific)
Molecular Weight:	371.10	Expiration/Reassay Date:	(batch specific)
CAS #:	4264-83-9 (Active Ingredient)	Storage:	
Appearance:		Grade:	REAGENT GRADE
Light yellow, smooth, flat tablets w/beveled edges		Additional Product Description:	
		Each 175 mg tablet contains 5 mg of p-Nitrophenyl Phosphate, Disodium. Balance is inert filler.	

---

**Additional Information**

TEST	SPECIFICATION	DISPOSITION
Abs.@309nm (1 tab/250mL Water)	0.48 - 0.58	PASS
Background Color	<= 0.14	PASS
Mean Tablet Weight	170 - 180 mg	PASS

Spec Set: 0405RG

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: \_\_\_\_\_

Title:

Date Printed: 05/04/2016

Page 1 of 1