



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

Corporate Headquarters

AMRESCO, LLC

28600 Fountain Parkway
Solon, OH 44139 USA
TEL: 440-349-1199 FAX: 440-349-1182
WWW.AMRESCO-INC.COM
Email: info@amresco-inc.com

CERTIFICATE OF ANALYSIS

Hydrochloric Acid

Code: 0369
Batch Number: 0926C125

Chemical Formula: HCl
Molecular Weight: 36.46
CAS #: 7647-01-0
Appearance:
Liquid, free of suspended matter

Manufacture Date: 04/26/2016
Expiration Date: 04/26/2017
Storage: Room Temperature
Grade: NF GRADE
Additional Product Description:
12N

TEST	SPECIFICATION	ANALYSIS	DISPOSITION
Bromide or Iodide	PASS / Fail	PASS	PASS
Free Bromine or Chlorine	PASS / Fail	PASS	PASS
Heavy Metals	<= 0.0005 %	<0.000005 %	PASS
Identification	PASS / Fail	PASS	PASS
Purity	36.5 - 38.0 %	37.5 %	PASS
Residue on Ignition	<= 0.008 %	<0.001 %	PASS
Sulfates Test	PASS / FAIL	PASS	PASS
Sulfites Test	PASS / FAIL	PASS	PASS

Spec Set: 0369NF

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Residual Solvents: No Class 1, 2, 3 or other residual solvents are present

Internal ID #: 546

Signature: Leona Edwardson
Title: Quality Control Sr. Manager

Date Printed: 11/22/2016

Page 1 of 1