



ISO 9001 CERTIFIED  
ISO 13485 CERTIFIED

## AMRESKO LLC

28600 Fountain Parkway  
Solon, Ohio USA 44139  
440/349-1199 FAX: 440/349-1182  
www.amresco-inc.com  
Email: info@amresco-inc.com

### CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

#### Erythromycin

**Code:** 0219

Chemical Formula:	C37H67NO13	Manufacture Date:	(batch specific)
Molecular Weight:	733.94	Expiration/Reassay Date:	(batch specific)
CAS #:	114-07-8		
Appearance:		Storage:	
White crystalline powder		Grade:	USP GRADE

#### Additional Information

TEST	SPECIFICATION	DISPOSITION
Crystallinity	PASS	PASS
Identification	PASS	PASS
Limit of Related Substances	PASS	PASS
Moisture (KF)	<= 10.0 %	PASS
Purity (Anhydrous)	85.0 - 100.5 %	PASS
Residue after Ignition	<= 0.2 %	PASS
Specific Rotation	-78 - -71 degrees	PASS
Thiocyanate	<= 0.3 %	PASS

Spec Set: 0219USP

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Residual Solvents, if applicable will be reported for actual batch number only.

Internal ID #: 269

Signature: \_\_\_\_\_

Title:

Date Printed:

02/29/2016

Page 1 of 1