



AMRESKO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Nalidixic Acid, Free Acid

Code: 0677

Chemical Formula:	C12H12N2O3	Manufacture Date:	(batch specific)
Molecular Weight:	232.24	Expiration/Reassay Date:	(batch specific)
CAS #:	389-08-2		
Appearance:		Storage:	
White to off-white crystalline powder		Grade:	HIGH PURITY GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Heavy Metals (as Pb)	<= 0.002 %	PASS
Loss on Drying	<= 0.5 %	PASS
Melting Point	225 - 231 Degrees C	PASS
Purity	>= 98 %	PASS
Residue on Ignition	<= 0.2 %	PASS

Spec Set: 0677HPG

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 269

Signature: _____
Title:

Date Printed: 05/13/2016
Page 1 of 1