



AMRESCO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Streptomycin Sulfate

Code: 0382

| | | | |
|-------------------|---|--------------------------|------------------|
| Chemical Formula: | (C ₂₁ H ₃₉ N ₇ O ₁₂) ₂ ·3H ₂ SO ₄ | Manufacture Date: | (batch specific) |
| Molecular Weight: | 1,457.38 | Expiration/Reassay Date: | (batch specific) |
| CAS #: | 3810-74-0 | | |
| Appearance: | | Storage: | |
| Fine white powder | | Grade: | USP GRADE |

Additional Information

| TEST | SPECIFICATION | DISPOSITION |
|----------------------|------------------|-------------|
| Identification | PASS | PASS |
| Loss on Drying (USP) | <= 5.0 % | PASS |
| pH (20%, Water) @25C | 4.5 - 7.0 | PASS |
| Potency (as is) | | PASS |
| Potency (Dry Basis) | 650 - 850 mcg/mg | PASS |

Spec Set: 0382USP

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Residual Solvents, if applicable will be reported for actual batch number only.

Internal ID #: 269

Signature: _____
Title:

Date Printed: 05/04/2016
Page 1 of 1