



AMRESCO LLC

28600 Fountain Parkway
Solon, Ohio USA 44139
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

CERTIFICATE OF QUALITY / CERTIFICATE OF ANALYSIS

Vancomycin Hydrochloride

Code: 0990

Chemical Formula:	C66H76Cl3N9O24	Manufacture Date:	(batch specific)
Molecular Weight:	1,485.67	Expiration/Reassay Date:	(batch specific)
CAS #:	1404-93-9		
Appearance:	White powder	Storage:	
		Grade:	USP GRADE

Additional Information

TEST	SPECIFICATION	DISPOSITION
Composition of Vancomycin: any other peak	<= 5.0 %	PASS
Composition of Vancomycin: vancomycin B	>= 85.0 %	PASS
Heavy Metals	<= 30 ppm	PASS
Identification	PASS	PASS
Moisture (KF)	<= 5.0 %	PASS
Monodechlorovancomycin	<= 4.7 %	PASS
pH (5%, Water) @25C	2.5 - 4.5	PASS
Potency (as is)		PASS
Potency (Dry Basis)	>= 900 mcg/mg	PASS

Spec Set: 0990USP

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

Residual Solvents, if applicable will be reported for actual batch number only.

Internal ID #: 269

Signature: _____

Title:

Date Printed:

06/09/2016

Page 1 of 1