


# 告別AMRESKO，正式進入VWR SYSTEM

## 波仕特感恩回饋四連發!



### (1) ACRYLAMIDE SOLUTION 買就送!

Cat No.	Product Name	特價
E344-500ML	ACRYL/BIS™ SOLUTION (30%) 29:1	\$2,280
E347-500ML	ACRYL/BIS™ SOLUTION (30%) 37.5:1	\$1,860

### 買就送TG BUFFER, 10X, 2LITER

### (2) VWR TRIS 震撼價!!

Tris, free base (#0497)  
5Kg \$10,200  
25Kg \$45,000



### (3) 高CP值 AGAROSE

VWR AGAROSE RA (#N605-100G)

特色: 高解析度, 低EEO, 高通量實驗最適用的Agarose!

### 超低價 NT\$900!!

### (4) PROZYME™ II DNA POLYMERASE

特色:

NT\$1360, 買一送一!



#PT-GL-TBR  
(250 units)

- 缺乏 3' → 5' exonuclease 活性
- 錯誤率不到 Taq DNA Polymerase 的二分之一
- 具 5' → 3' exonuclease 活性, 可用於 nick-translation 實驗
- 由Thermus Brockianus (Tbr) 純化出來, 熱穩定性高
- 含10X Reaction Buffer: 100 mM Tris-HCl · pH 8.8, 500 mM KCl · 15 mM MgCl<sub>2</sub> · 1% Triton X-100

另有QR  
神秘價

