


新年新春新品來報 實驗順利揚眉兔氣

6X loading dye 安全核酸染劑

藍光或 UV 照膠儀器皆適用，可選擇符合習慣、最合適的品項！

SafeDye for Blue Light (6X)

- #PT-D3002, 1ML
- 靈敏度: 0.3 ng of DNA per gel band
- 使用藍光，降低DNA的斷裂機率
- 偵測 dsDNA, ssDNA and RNA in agarose and polyacrylamide gels




PT-D3002 競爭品 PT-D3002 競爭品

DNA ladder (#PT-M5-1KB) were electrophoresed through 1% agarose gel (#HI-BIOMAX-100) with 1-fold and 10-fold dilution. Each lane of 5 μ L ladder were stain with 1 μ L SafeDye for Blue Light (6X) (#PT-D3002).

SafeDye for UV (6X)

- #PT-D4001, 1ML
- 靈敏度: 0.1 ng of DNA per gel band
- DNA 片段飄移的現象較輕微
- 偵測 dsDNA, ssDNA and RNA in agarose and polyacrylamide gels



PT-D4001 競爭品 PT-D4001 競爭品

DNA ladder (#PT-M5-1KB) were electrophoresed through 1% agarose gel (#HI-BIOMAX-100) with 1-fold and 10-fold dilution. Each lane of 5 μ L ladder were stain with 1 μ L SafeDye for UV (6X) (#PT-D4001).

買安全染劑
送Marker

新春限定

新品買一送一
或好禮五選一

TAE 1L

TBE 1L

電影票 1 張

50bp ladder #PT-M1-50

100bp ladder #PT-M1-100T

PT-M1-50 PT-M1-100T



0800-231-530

加碼特惠

20% OFF

兩組更優惠
- 快來洽詢 -

Gene-Spin™ PCR Clean up & Gel Extraction Kit

#Cat. PT-CG250

- Spin Column 方式純化
- 純化 100 bp ~ 10 kbp 的 DNA 片段
- 回收率: PCR clean-up > 80%,
Gel Extraction > 60 %
- Binding capacity: ≤ 25 µg DNA

\\ 操作時間最快只要20分鐘 /

