

Acry/Bis-Acrylamide 29:1 30% (w/v) Solution

Certificate of Analysis

Cat. No. : PT-BF344-500ML

Lot No. :

Reassay Date:

Quantity: 500 mL

Grade: UltraPure

Analysis:

Item	Specification	Results
Appearance	Clear, colorless solution	Pass
Concentration of acrylamide (w/v)	29.0 %	Pass
Concentration of bisacrylamide (w/v)	1.0 %	Pass
Conductivity (@18°C)	≤ 5 μmho	1.61 μmho
DNase, RNase and Protease	None detected	Pass
Heavy metals	<0.005%	Pass
Fe, Pb	<0.001%	Pass
Conclusion	This results conforms with enterprise standards	

Approved by: _____

Title: QC SUPERVISOR

Date: _____

For Research Using Only

Manufactured for and distributed by Protech Technology Enterprise Co., Ltd.

TEL: +886-2-2655-7677 / FAX: +886-2-2655-7601 / Toll Free: 0800-231-530

E-MAIL: service@bio-protech.com.tw; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>