
eXTReme™ SealMate™ Rolls

Ready to roll with eXTReme efficiency? eXTReme™ Seal features an extra-strong thick adhesive to minimize evaporation yet removes easily without leaving residue on the plate.

Application:
PCR

Temperature Range:
-20°C to +120°C




CATALOG #	DESCRIPTION	QUANTITY
SM-KIT-XTR	eXTReme Seal Starter Kit Dispenser + 2 Rolls	80/Roll
RSMXTR-2	eXTReme Seal Film 2 Rolls (Refill)	80/Roll


Features and Benefits:

- Certified DNase- & RNase- free
- Extra-strong thick adhesive minimizes evaporation, yet removes easily without leaving residue on the plate
- Environmentally friendly; no paper liner to remove and discard
- Each sheet is separated by an adhesive-free, perforated, color-coded blue band
- The smooth lower surface of the dispenser can be used to seal the film to the plate
- Dispenser protects unused seals from contamination
- Lower support prop prevents dispensing bill from contacting surfaces to prevent contamination
- Replacement rolls can be easily installed by releasing the locking clasps on the dispenser
- Compatible with flat-top PCR plates

The SealMate™ Dispenser Makes Microplate Sealing Easy!



Open the dispenser to insert a new roll



Close and lock the side latches to dispense

Excel Scientific
18350 George Blvd.
Victorville, CA 92394 USA

t +1.760.246.4545
info@excelscientific.com
excelscientific.com

EXCEL SCIENTIFIC